Veroval®

2 in 1 ECG and Blood Pressure Monitor

**FOR STROKE
PREVENTION**

The Veroval® ECG and Blood Pressure Monitor combined device is a single channel ECG that can detect various arrhythmias, for example atrial fibrillation, extrasystoles and therefore helps to make a preliminary assessment by carrying out a basic recording.

OK

Pocket ECG and Blood Pressure Monitor to record the problem the moment it happens

**For stroke
prevention**

Portable ECG function:

- ✓ Detects atrial fibrillation – a major cause of strokes
- ✓ Detects cardiac arrhythmias, such as extrasystoles – an early sign of circulatory disorders and inflammation of the heart muscle
- ✓ Recording of the heart rhythm in 30 seconds – at anytime and anywhere
- ✓ Immediate, easy to understand results display as result indicators or as a graph for doctors

Fully automatic blood pressure measurement:

- ✓ High measuring accuracy
- ✓ Gentle and precise measurement on the upper arm with a comfortable and ergonomic cuff
- ✓ Simple colour-coded classification of results
- ✓ Memory for 128 readings (64 readings per user, max. 2 users)
- ✓ Developed in Germany

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TRACK YOUR HEART
with **Veroval®** ECG and
blood pressure monitor

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Did you know that
Stroke is #2 cause of
death in Oman, #3 in
Kuwait and Qatar?¹

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Atrial fibrillation:

You may have it and not know it!²

It's a type of arrhythmia that can be
silent or have a **triad of
symptoms**.³

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


Why atrial fibrillation is dangerous?

It can lead to blood clots causing stroke and can only
be detected when it happens!⁸

It can be asymptomatic or characterized by a markedly
accelerated, constantly irregular heartbeat.⁹

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
50+

Having or not Hypertension Cardiac Arrhythmias

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ECG recording as per your doctor's advice, for silent AFib, or for heart palpitations, restlessness, perspiration and breathing difficulties.

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What is ECG?

A contraction of the heart muscle is caused by an electrical stimulus that is detected on the surface of the body and displayed by an ECG graphic.

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تذكر! Keep in mind!

As a customer, do not interpret the electrocardiograms by yourself and do not use them for self-prescribed-treatment.

Recording a 1-lead ECG is not a substitute for medical treatment!

The device is not a substitute for either a medical examination of the heart function, nor recording of a medical ECG that must be obtained using more elaborate measuring equipment.

The device does not diagnose any medical conditions that may be the possible cause of the change observed in the ECG.

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
Method 1: Chest

Has a high measurement stability.

Hold the device in your right hand. Place your right index finger on the upper electrodes.

Place the lower electrode on bare skin about 5 cm below your left nipple.

Ensure that there is no skin contact between your hand and chest or the recording will not be correct.

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
Ensure that there is no skin contact between your right and left hands or the recording will not be correct.

Press with moderation and steadiness during the recording.

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Place your right index finger on the upper electrodes.

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
The ECG recording can be invalid for the following reasons:

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If the heart rate flashes, this means that the ECG signals are unstable or weak. In this case, repeat the recording.

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If you cannot record in method 1, or the ECG is frequently showing unstable **EE** symbol, change to method 2.

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

If the contact between the electrodes and the finger is not sufficiently firm.

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
If you cannot record any stable measurements, it may simply be due to harmless causes such as the heart shape or due to pathological reasons which needs a doctor advice.

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Sources:

World Health Data Platform „Global Health Estimates: Life expectancy and leading causes of death and disability - Kuwait 2019“, accessed: 2021-01-28, www.who.int/data/gho/data/themes/mortality-and-global-health-estimates/gh-leading-causes-of-death, WHO World Health Data Platform „Global Health Estimates: Life expectancy and leading causes of death and disability - Oman 2019“, accessed: 2021-01-28, www.who.int/data/gho/data/themes/mortality-and-global-health-estimates/gh-leading-causes-of-death. WHO World Health Data Platform „Global Health Estimates: Life expectancy and leading causes of death and disability - Qatar 2019“, accessed: 2021-02-09, www.who.int/data/gho/data/themes/mortality-and-global-health-estimates/gh-leading-causes-of-death. 2. Sunrise Hospital & Medical Center (2017-01-23) „Atrial fibrillation: You may have it and not know it“, accessed: 2021-01-28, www.lasvegassun.com/hative/sunrise/2017/jan/23/atrial-fibrillation-you-may-have-it-and-not-know-it/. 3. American Heart Association (c2021) „FAQs of Atrial Fibrillation“ accessed: 2021-01-28, www.heart.org/en/health-topics/atrial-fibrillation/afib-resources-for-patients-professionals/afib-awareness. 4. Rienstra, Michiel et al. „Symptoms and functional status of patients with atrial fibrillation: state of the art and future research opportunities.“ Circulation vol. 125,23 (2012): 2933-43. doi:10.1161/CIRCULATIONAHA.111.069450. 5. Boriani, Giuseppe et al. „Asymptomatic atrial fibrillation: clinical correlates, management, and outcomes in the EORP-AF Pilot General Registry.“ The American journal of medicine vol. 128,5 (2015): 509-18.e2. doi:10.1016/j.amjmed.2014.11.026. 6. Go, A S et al. „Prevalence of diagnosed atrial fibrillation in adults: national implications for rhythm management and stroke prevention: the Anticoagulation and Risk Factors in Atrial Fibrillation (ATRIA) Study.“ JAMA vol. 285,18 (2001): 2370-5. doi:10.1001/jama.285.18.2370. 7. Miyasaka, Yoko et al. „Secular trends in incidence of atrial fibrillation in Olmsted County, Minnesota, 1980 to 2000, and implications on the projections for future prevalence.“ Circulation vol. 114,2 (2006): 119-25. doi:10.1161/CIRCULATIONAHA.105.595140. 8. ARNI Institute, T. Balchin (2018) „Atrial Fibrillation raises risk of stroke by 5: Can you tell if you have it?“, accessed: 2021-01-28, www.arni.uk.com/atrial-fibrillation-raises-risk-of-stroke-by-5-can-you-tell-if-you-have-it/. 9. American Heart Association (c2019). „AFib Can Happen to Anyone“, accessed: 2021-01-28, www.heart.org/en/health-topics/atrial-fibrillation/afib-resources-for-patients-professionals/afib-can-happen-to-anyone-infographic.

PAUL HARTMANN Middle East
Dubai Airport Freezone
P.O. Box 54525 Dubai
United Arab Emirates
Phone: +971 4 2996996

me.hartmann.info