5.1 Upper arm blood pressure monitor


Dear customer
 Thank you for choosing a HARTMANN blood pressure monitor. The Veroyal® BPU 22 is a quality product intended for fully automatic blood pressure measurement on the upper arm of adults and is suitable for clinical and domestic use. Requiring no preliminary settings, this blood pressure monitor easily and automatically inflates for simple, quick and reliable measurement of the systolic and diastolic blood pressure and the pulse rate.
 The monitor also indicates if it detects an irregular heartbeat.

We wish you the best of health.
 Please read these instructions for use carefully before handling the device for the first time. Correctly measuring your blood pressure depends on an appropriate handling of the device. These instructions for use are designed to teach you each step in taking your own blood pressure measurements using the Veroyal® BPU 22. You will be given important and useful hints to ensure that you obtain reliable results for your personal blood pressure profile. Use the monitor in accordance with the information provided in the instructions for use. Please store these instructions carefully and ensure other users can access them. Check the packaging for any damage and that the contents are complete.


- Package contents**
- Blood pressure monitor
 - Universal upper arm cuff
 - 4 x 1.5 V AA batteries
 - Storage bag
 - Instructions for use with warranty certificate

1. Device and display description


- Blood pressure monitor**
- 1 Extra large LCD display
 - 2 START/STOP button
 - 3 Memory button for User 1
 - 4 Memory button for User 2
 - 5 Cuff connector socket
 - 6 Battery compartment

- Cuff**
- 7 Cuff with application instructions
 - 8 Size indicator to correctly fit cuff
 - 9 Cuff plug
 - 10 Cuff tube

- Display**
- 11 Systolic blood pressure
 - 12 Diastolic blood pressure
 - 13 Pulse rate
 - 14 Flashes when the device is measuring and the pulse rate is detected
 - 15 Irregular heartbeat symbol
 - 16 Cuff fit check symbol
 - 17 Battery symbol
 - 18 Traffic light system for your values/progress bar
 - 19 Average value (AM), morning (AM), evening (PM) / memory space number
 - 20 User memory
 - 21 Date and time display

2. Important information

See table on extra page.

Important instructions for use

Intended use
 The fully automatic blood pressure monitor Veroyal® BPU 22 is intended for clinical and home use for the non-invasive measurement of diastolic and systolic blood pressures and pulse rate of adults by means of an inflatable cuff which is wrapped around the upper arm.
 The use of the device does not require a specific knowledge or professional ability.
Intended user
 The patient is the intended operator, except in case of patients that require special assistance.
Patient population
 • Adults whose upper arm circumference is within the range indicated on the cuff.
 • Only use the monitor to take blood pressure measurements on the upper arm. Do not place the cuff on other areas of the body.
 • Only use the supplied cuff or an original replacement cuff. The measured values may otherwise be incorrect.
 • If measured values are questionable, repeat the measurement.

- Do not leave the monitor unattended near small children or persons who cannot operate it themselves. There is a risk of strangulation if the cuff tube becomes wrapped around the neck. Swallowing small parts that have detached from the device can also result in suffocation.
- Do not under any circumstances carry out blood pressure measurements on newborns, babies or toddlers.
- Do not use the monitor for diagnosis in clinical situations where the patient is in immediate danger.
- Please note that the pressure build-up in the cuff can lead to temporary disruption of medical devices being used at the same time on the same arm.
- The function of the arm on which the cuff is placed may become impaired during inflation.
- If you are carrying out the measurement on another person, please ensure that using the blood pressure monitor does not have a lasting negative effect on the blood circulation.
- People with restricted physical, sensory or mental skills should be supervised by a person responsible for their safety and receive instructions from this person on how to use the device.
- Taking measurements too often within a short time period as well as maintaining continuous cuff pressure can disrupt the blood circulation and cause injuries. Please take a break between measurements and do not bend or fold the air tube. If the monitor malfunctions, remove the cuff from the arm.
- Always consult your doctor about if and when the blood pressure monitor can be used by pre-eclampsia patients during pregnancy.
- Always consult your doctor about if and when the blood pressure monitor can be used by persons who have had a mastectomy or lymph node clearance.
- If a possible allergic skin reaction occurs on the arm where the cuff is used, discontinue its use and consult a physician.
- Do not use the blood pressure monitor together with a high frequency surgical device.
- The device is not intended for use in vehicles (e.g. ambulances) or helicopters.

- Contraindications**
- Please do not place the cuff over a wound as this may cause further injuries.
 - If an intravenous treatment is being carried out or a venous catheter is inserted in the arm, blood pressure measurements can lead to injury. Never use the cuff on the arm to which these conditions apply.
 - If you are prone to hematomas, do not put on a cuff nor carry out a blood pressure measurement.
 - Do not use the blood pressure monitor in the absence of a full arm.

- Important instructions for self-measurements**
- Even minor changes in internal and external factors (e.g. deep breathing, stimulants, talking, excitement, environmental factors) can lead to fluctuations in blood pressure. This explains why deviating values are often measured at the doctor or pharmacist.
 - Measurement results fundamentally depend on the measurement location and the position of the patient (sitting, standing, reclining). They are also influenced by e.g. physical activity and physiological preconditions of the patient. For comparable values, carry out the measurement in the same location and in the same position.
 - Cardiovascular disease can result in erroneous readings or lower the accuracy of the measurement. The accuracy may also be affected if you have very low blood pressure, diabetes, circulatory disorders and arrhythmias and chills or tremors.

- Consult your doctor before measuring your own blood pressure if ...**
- you are pregnant. Blood pressure may vary during pregnancy. Regular blood pressure monitoring is particularly important if you have high blood pressure, because elevated values may affect foetal development. Always check with your doctor if and when you should measure your own blood pressure, especially if you suffer from pre-eclampsia.
 - you have diabetes, liver function disorders or narrowing of blood vessels (e.g. arteriosclerosis, peripheral arterial disease) as these may cause anomalous values.
 - you have certain blood disorders (e.g. haemophilia) or severe circulatory disorders or are taking blood-thinning medication.
 - you are on dialysis or taking anticoagulant medication (anticoagulants), platelet aggregation inhibitors or steroids.
 - you wear a pacemaker, which may cause anomalous values. The blood pressure monitor itself does not affect a pacemaker. Please note that the pulse rate displayed is not suitable for checking the frequency of pacemakers.
 - you tend to bruise easily and/or are sensitive to pressure pain.
 - you suffer from severe heart rhythm disorders or cardiac arrhythmias. Due to the oscillometric measurement method used, incorrect measured values may be obtained in some cases or it may not be possible to measure blood pressure.
 - this symbol appears often, it may indicate heart rhythm disorders. Please consult your doctor in that case. Severe heart rhythm disorders may sometimes produce false measurements or reduce the measurement accuracy. Please discuss with your doctor whether taking your own blood pressure is appropriate for you.

- The readings taken by you are for your information only – they are not a substitute for a medical examination. Discuss your measured values with your doctor and never make a medical decision based on these measurements (e.g. medications and their dosage).
- Taking your own blood pressure is not a substitute for medical treatment. Do not interpret your measured values yourself and do not use them for self-prescribing treatment or for self-prescribing your own medication. Do not trust his or her diagnosis. Take medication only as prescribed by your doctor and do not alter the dose on your own accord. Discuss with your doctor the appropriate time to measure your blood pressure.

An irregular heartbeat is present when the heart rhythm deviates by more than 25 % from the median heart rhythm. The contraction of the heart muscle is stimulated by electrical signals. A disorder of these electrical signals is called an arrhythmia. This may be due to a genetic predisposition, stress, age, lack of sleep, exhaustion or similar factors. A doctor can establish whether an irregular heartbeat is the result of an arrhythmia.

Heart rhythm disorders are disturbances of the normal rate or rhythm of the heartbeat. A distinction should be made between mild and severe heart rhythm disorders. This can only be determined by a special examination carried out by a doctor.

- Power supply (batteries, power supply unit)**
- Observe the polarity labels plus (+) and minus (-).
 - Only use high-quality batteries (see specification in Chapter 11 'Technical data'). If you use low-quality batteries, we can no longer guarantee the specified number of measurements.
 - Never mix old and new batteries or batteries from different manufacturers.
 - Remove flat batteries immediately.
 - Replace batteries if the battery symbol remains illuminated.
 - Always replace all the batteries at the same time.
 - If you do not intend to use the monitor for a longer period of time, you should remove the batteries to prevent possible leakage.
 - If you use the monitor with the power supply unit, position the monitor so that you can disconnect the power supply at any time.

- Information about batteries**
- **Choking hazard**
 Small children could swallow batteries and suffocate on them. Keep the packaging, batteries and the monitor out of the reach of children.
 - **Risk of explosion**
 Do not throw batteries into a fire.
 Batteries must not be charged or short-circuited.
 If a battery has leaked, wear protective gloves and clean the battery compartment with a dry cloth. If liquid from a battery cell comes into contact with skin or eyes, clean the affected area with water and seek medical attention if necessary.
 - Protect batteries from excessive heat.
 - Do not disassemble, open or crush batteries.

- Safety information concerning the device**
- This blood pressure monitor is not waterproof.
 - This blood pressure monitor is made of high-quality electronic precision components. The accuracy of the measured values and the lifetime of the device depend on careful handling.
 - Protect the device from violent shaking, hitting or vibrations and do not let it drop on the ground.
 - Do not bend or fold the cuff excessively.
 - Never open the device. The device must not be modified or dismantled, and must not be repaired by the user. Repairs may only be carried out by authorized specialists.
 - Do not inflate the cuff when it is not properly applied to the upper arm.
 - Use the device only with the approved upper arm cuff. The device may otherwise sustain external or internal damage.

- The cuff tube may only be removed from the unit by pulling on the corresponding connector plug. Never pull on the tube itself.
- Do not expose the device to extreme temperatures, humidity, dust or direct sunlight as this may cause it to malfunction.
- Keep the packaging, batteries and device out of reach of children.
- Protect the device and cuff against contact with pets and pests in order to avoid damage.
- Please comply with the storage and operating conditions defined in Chapter 11 'Technical data'. Storing or using the device outside the specified temperature and humidity range can affect the accuracy of the measurement or the function of the device.
- If the device was not stored within the minimum/maximum permissible storage conditions, a waiting period of at least 2 hours must be observed before using it under the specified operating conditions (Chapter 11 'Technical data') or an ambient temperature of approx. 20 °C.
- Do not use the device in explosive environments with inflammable gases or concentrated oxygen.

- Instructions regarding electromagnetic compatibility**
- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
 - It may only be possible to use the device to a limited extent in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.
 - Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. However, if it is necessary to use the device in the manner described above, this device and the other devices must be monitored to ensure that they are working properly.
 - The use of accessories other than those specified or provided by the manufacturer of this device may result in increased electromagnetic emissions or a decrease in the device's electromagnetic immunity and lead to faulty operation.
 - Failure to comply with the above can impair the performance characteristics of the device.
 - Do not use the device near strong electromagnetic fields and keep it away from radio equipment and mobile phones. Portable and mobile high-frequency and communication devices, such as telephones and mobile phones, can impair the functionality of this electronic medical device. Therefore keep a minimum distance of 30 cm.
 - Further documentation in accordance with the electromagnetic compatibility standard can be requested from the manufacturer or the customer service (contact details see Chapter 10).

Instructions for the measurement function check
 Every Veroyal® device has been carefully tested by HARTMANN for measurement accuracy and has been developed with a view to a long service life. We recommend carrying out a metrological check every 2 years for devices in professional use, for example, in pharmacies, medical practices or hospitals. Please also observe the national regulations specified by the legislator. Metrological checks may only be carried out by the competent authorities or authorized maintenance providers against cost reimbursement.

- Disposal information**
- This product (monitor and cuff) is subject to the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment and is marked accordingly. Never dispose of electronic equipment in your household waste. Please obtain information about your local regulations on the proper disposal. This protects the environment and human health.
 - The batteries comply with the European Directive 2013/56/EU. To protect our environment, empty batteries must not be disposed of in household waste. Please comply with the relevant waste-disposal regulations or use public collection points, recycling points or electronic retailers.

3. Blood pressure information

- Two values must be measured to determine blood pressure:
- Systolic (upper) blood pressure, which is produced when the heart contracts and pumps blood into the blood vessels.
- Diastolic (lower) blood pressure, which is produced when the heart is expanded and fills with blood again.
- Blood pressure values are displayed in mmHg.

For improved assessment of the results, a traffic light system is located on the left side of the Veroyal® BPU 22 which is a direct result indicator, making it easier to categorise the measured value. The World Health Organisation (WHO) and the European Society of Cardiology (ESC) & European Society of Hypertension (ESH) have developed the following summary for classifying blood pressure values:

Resulting indicator	Evaluation	Systolic pressure	and/or	Diastolic pressure
Red	Grade 3 hypertension	above 179 mmHg	and/or	above 109 mmHg
Orange	Grade 2 hypertension	160 – 179 mmHg	and/or	100 – 109 mmHg
Yellow	Grade 1 hypertension	140 – 159 mmHg	and/or	90 – 99 mmHg
Green	High normal	130 – 139 mmHg	and/or	85 – 89 mmHg
Green	Normal	120 – 129 mmHg	and/or	80 – 84 mmHg
Green	Optimal	up to 119 mmHg	and	up to 79mmHg

- Established hypertension (high blood pressure) is defined as having a systolic value of at least 140 mmHg and/or a diastolic value of at least 90 mmHg.
- In general, low blood pressure (hypotension) is defined as a systolic value below 105 mmHg and a diastolic value of 60 mmHg. However, the threshold between normal and low blood pressure (hypotension) is not as accurately specified as the threshold for high blood pressure (hypertension). Hypotension may be associated with symptoms such as dizziness, fatigue, tendency to faint, visual disturbances or high pulse rate. In order to ensure that hypotension or the corresponding symptoms are not attendant symptoms of serious illness, a doctor should be consulted in case of doubt.


Constantly elevated blood pressure significantly increases the risk of developing other disorders. Physical consequences of high blood pressure, e.g. heart attack, stroke and organ damage, are the most common causes of death worldwide. Daily blood pressure monitoring is therefore an important measure that can help to protect you from these risks. You should discuss this with your doctor, especially if you have frequently elevated or borderline blood pressure values (see Chapter 3 'Blood pressure information'). Your doctor will initiate appropriate measures.

4. Preparing the measurement


- Inserting/changing the batteries**
- Open the battery cover on the bottom of the device. Insert the batteries (see Chapter 11 'Technical data'). Ensure correct polarity (+ and -) when inserting the batteries. Close the battery lid. or flash on the display. Set the date and time as described above.
 - If the Change battery symbol is permanently illuminated, blood pressure can no longer be measured and you need to replace all the batteries.

Setting the date, time
 It is very important that you set the date and time correctly. This is the only way to save your measured values correctly, with the date and time, for subsequent retrieval. This is also important for correctly using all the memory and analysis functions.

- To switch to settings, reinsert the batteries or press the START/STOP button for 5 seconds. Proceed as follows:

- If the 12-hour format is selected, the month is displayed before the day.

The device switches off automatically once all data has been set.


The measured values will remain in the memory when the batteries are changed. The date and time settings will also be saved.

5. Measuring your blood pressure

The 10 golden rules of blood pressure measurement
 Many factors play a role in blood pressure measurement. These ten general rules help you to correctly carry out the measurement.

1. Rest for about 5 minutes before measuring your blood pressure. Even working at your desk increases blood pressure by an average of about 6 mmHg systolic and 5 mmHg diastolic.
2. Do not consume any nicotine or caffeine up to one hour before measurement.
3. Do not measure when you have a strong urge to urinate. A full bladder can increase blood pressure by about 10 mmHg.
4. Measure your blood pressure on a completely bare upper arm and in a comfortable upright seated position. The blood circulation must not be constricted by a rolled-up sleeve, for example.
5. If you are using a wrist monitor, hold the cuff at heart height during measurement. The cuff of an upper arm monitor is naturally positioned at the correct height.
6. Please do not talk or move when measuring your blood pressure. Talking increases values by about 6–7 mmHg.
7. Wait for at least one minute between each measurement so that the blood vessels are ready, by removing the pressure applied to them.
8. The date and time are conveniently and easily recorded with the measured values as well as any medications you have taken, with the Veroyal® medic.connect app.
9. Take measurements regularly. Even if your values have improved, you should continue to check them for monitoring purposes.
10. Always take measurements at the same time of day. Because a person has about 100,000 different blood pressure values every day, individual measurements are not meaningful. Only regular measurements taken at the same time each day over a longer period of time give a meaningful assessment of blood pressure values.

- Applying the blood pressure monitor**
- Before applying the cuff, insert the cuff socket into the connector socket on the left side of the device.
 - Do not mechanically constrict, compress or bend the sleeve tube.
 - Blood pressure must be measured on a bare upper arm. If the cuff is completely open, thread the end of the cuff through the metal ring to form a loop. Position the cuff on the upper arm so that the lower edge is 2–3 cm above the crook of the arm and above the artery. The tube points to the centre of the palm.

- Now bend your arm slightly, grasp the free end of the cuff, wrap it firmly round your arm and close the hook-and-loop fastener.

The cuff should fit firmly but not tightly. You should be able to insert two fingers between the arm and the cuff. Make sure the tube is not twisted or damaged.

Important: It is essential that the cuff is correctly applied to obtain correct measurements. The cuff included with the device is a universal cuff for upper arm circumferences of 22 to 42 cm. The white arrow must point to an area within the sizing scale. If it is outside the range, it is not possible to guarantee a correct measurement result.

This innovative Veroyal® device with Comfort Air Technology ensures comfortable measurement. The individual inflation pressure is determined independently for each measurement and depends on the corresponding systolic blood pressure value. This allows for more comfortable upper arm measurements.


- Carrying out the measurement**
- You should take your blood pressure in a quiet place, in a relaxed and comfortable seated position.
 - Measurement can be taken on the right or left arm. We recommend you perform the measurement on the left upper arm. Over the long term, the arm giving higher readings should be used for blood pressure monitoring. However, if there is a very clear difference between readings on either arm, you need to check with your doctor which arm you should use for the measurement.
 - Always measure on the same arm and place your forearm in a relaxed position on a support.
 - We recommend that you sit with your back against the back of the chair. Place both feet flat on the floor next to each other. Do not cross your legs. Relax your forearm and hand with your palm facing upwards on a support and make sure the cuff is at the height of your heart.
 - Do not measure your blood pressure after taking a bath or playing sports.

- Do not eat, drink or exercise for at least 30 minutes before the measurement.
- Please wait for at least one minute between two measurements.

Starting the measurement

- Start a measurement only after applying the cuff because the cuff may otherwise be damaged by the excess pressure that is created. Press the START/STOP button . The appearance of all display segments, followed by the time and date, indicates that the device is carrying out an automated check and is ready for use.
- Check the display segments for completion (see Chapter 1 'Device and display description').
- After about 0.5 seconds the cuff inflates automatically. If this inflation pressure is insufficient or if interference is experienced during measurement, the device will increase pressure to the appropriate higher pressure value in increments of 40 mmHg. During inflation, the results indicator on the left of the display also increases at the same time.
- If the cuff has been applied firmly enough on the arm, the cuff symbol appears on the display. If the cuff symbol does not appear in the display, the cuff has not been fitted securely enough and the error message appears on the display after a few seconds.

- Do not speak or move for the duration of the measurement.
- As the cuff pressure is released, the heart symbol flashes and the dropping pressure in the cuff is shown on the display.
- After the measurement is complete, the systolic and diastolic blood pressure as well as the pulse rate appear simultaneously.

- In addition to the measured values, the time, date, the associated user memory or as well as the associated memory numbers (e.g. M05) are shown. While the result of the measurement is displayed, you can assign the values to the corresponding user memory by pressing the button. If you do not assign the values to the corresponding user memory, the measured values are automatically stored in the user memory shown on the display.
- You can use this results indicator on the left of the display to classify your measurement result (see table in Chapter 3 'Blood pressure information').
- To switch off the device, press the START/STOP button , otherwise the device will automatically switch off after 1 minute.

- If you wish to stop the measurement for any reason, simply press the START/STOP button . The inflation or measurement process is interrupted and cuff pressure is automatically released.
- If an irregular heartbeat was detected during the measurement, the display will show the heart symbol . However, the measurement may also have been disrupted by body movement or speaking. It is best to repeat the measurement. If you see this symbol regularly when measuring your blood pressure, we recommend that your doctor checks your heart rhythm.

6. Memory function

- User memory**
- The Veroyal® BPU 22 stores up to 100 measurements in each user memory. The most recent measurement is always saved together with the date and time in memory slot no. 1 and all older saved values are moved up one memory slot. Once all memory slots are occupied, the oldest value is deleted.
 - The device has two memory buttons, and , that are used to record the measurement results for two separate users. represents the measured values for the first user, represents the measured values for the second user. After the measurement is finished, you can assign the measured value to the respective person by pressing and . This assignment can be made for 10 seconds. If you do not assign them, the measured value is automatically stored in the user memory that is currently displayed.
 - The time of the measurement is also saved, along with each of the blood pressure values, to enable the morning and evening averages to be determined. The time saved in the device must therefore be set correctly to the actual time (see Chapter 4, section 'Setting the date and time').


Veroyal® BPU 22 features the following measurement data storage functions (according to the guidelines of the European Society of Hypertension (ESH)):

- Individual measured value memory
- Average value of all measured blood pressure values per user
- Average value of the morning values
- Average value of the evening values



- The memory recall is activated by pressing the or button with the device switched off. For values saved in the first user memory button, press the button, for the second user memory press the memory button.

Average values

- After selecting the respective user memory, the corresponding symbol or and appear on the display. The average value of all the data stored for the particular user memory is displayed.

- By pressing the button again (or the button if you are in user memory 2), the average values for all morning measurements 'PMT' (00:00 to 11:59) for the past 7 days (including the current day) appear.
- By pressing the button again (or the button if you are in user memory 2) the average values for all evening measurements 'PNT' (12:00 to 23:59) taken on the last 7 days (including the current day) appear.

- If there is no morning or evening value for the past 7 days in the memory, dashes will be shown on the display instead of the morning or evening average values. If no values at all are stored in the memory, then dashes will also be displayed instead of the overall average value.

- Individual measured values**
- By pressing the button again (or the button if you are in user memory 2), all memory values can be retrieved one after the other, starting with the most recently measured value.
 - If an irregular heartbeat was detected during measurement, that information is also stored and displayed when the measured value is retrieved from the device's memory along with the systolic and diastolic blood pressure readings, pulse rate, time, date and year.
 - You can cancel the memory recall at any time by pressing the START/STOP button .
 - Otherwise, the device will switch off automatically after a few seconds.
 - Even if the power supply is interrupted, e.g. when changing the batteries, the stored values will still be available.

Deleting saved readings
 All data saved can be deleted separately for user memory and user memory . To do so, press the button of the corresponding user memory (or). The average 'A' value will then appear on the display. Press and hold the user memory button for 5 seconds. 'CL 00' appears on the display. All data in the selected user memory are now deleted. If you release the button early, no data will be deleted.

7. Explanation of error messages

Error	Possible causes	Remedy
Monitor cannot be switched on	Batteries are missing, inserted incorrectly or flat.	Check the batteries. If necessary, replace with four new identical batteries.
Cuff will not inflate	The cuff connector plug is incorrectly positioned in the connector socket of the monitor. Wrong cuff type connected.	Check the connection between the cuff connector plug and socket. Verify that only an approved Veroyal® cuff and connector are being used.
E1	Pulse could not be detected correctly.	Check that the cuff has been positioned correctly. Do not talk or move during measurement.
E2	Systolic or diastolic blood pressure could not be detected.	Check that the cuff has been positioned correctly. Do not talk or move during measurement.
E3	The cuff was applied too firmly or too loosely or the cuff is defective.	Apply the cuff so that two fingers fit between the cuff and the upper arm.
E4	The air tube is not correctly connected to the monitor. Check that the connector plug is correctly positioned. If this error occurs frequently, you should use a new cuff.	
E5	Error during measurement.	Please rest for at least 1 minute, then take the measurement again (see Chapter 5 'Measuring your blood pressure').
E6	Inflation pressure is higher than 300 mmHg.	Please rest for at least 1 minute, then take the measurement again (see Chapter 5 'Measuring your blood pressure').
E7	A system error has occurred.	Contact customer service if this error message appears.
E8	The batteries are almost empty.	Replace batteries.
Implausible measured values	Implausible measured values often occur due to inappropriate handling of the device or mistakes during the measuring procedure.	Please refer to Chapter 5 'Measuring your blood pressure' as well as the safety information. Then repeat the measurement.

- Switch the device off if an error symbol appears. Check for all possible causes and follow the 10 golden rules (Chapter 5 'Measuring your blood pressure') as well as the instructions for self-measurement from Chapter 2 'Important information'.
- Relax for one minute and then take the measurement again.

- 8. Monitor maintenance**
- Only use a soft, damp cloth to clean the device. Please do not use thinner, alcohol, detergents or solvents.
 - The cuff can be cleaned carefully using a slightly damp cloth and mild soap solution. Do not completely immerse the cuff in water.
 - It is recommended to clean and disinfect the cuff regularly or after each use, especially when used by several users, to prevent infection. The cuff should be disinfected, particularly on the inside, by wiping with a disinfectant. Use a disinfectant that is compatible with the cuff materials, e.g. 75 % ethanol or isopropyl alcohol. To protect from external influences, keep the device, the cuff and these instructions in the storage bag.
 - If you store a device and cuff, do not place any heavy objects on the device and the cuff. Remove the batteries.

- 9. Warranty conditions**
- We offer a 3-year warranty for this high-quality blood pressure monitor from the purchase date and in accordance with the conditions listed below.
 - Warranty claims must be made during the warranty period. The purchase date must be documented by an appropriately completed and stamped warranty certificate or proof of purchase.
 - Within the warranty period, HARTMANN will replace or repair free of charge any faulty device components caused by material faults or manufacturing errors. This does not extend the warranty period.
 - For local warranty conditions, please consult the local warranty certificate.
 - This device is solely intended for the purpose stated in these instructions for use.
 - The manufacturer is not liable for damage or injuries resulting from improper handling, use or unauthorized interference. Accessory parts that are subject to wear and tear (batteries, cuffs, etc.) are excluded from the warranty. Claims for compensation are limited to the value of the product; compensation for consequential damage or injury is expressly excluded.
 - In warranty cases, please send the monitor with cuff and, if applicable, the power

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