## How to use?

There are 2 methods to record the ECG:

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### Method 1: Chest

Has a high measurement stability.

Hold the device in your right hand. Place your right index finger on the upper electrodes.

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
Place the lower electrode on bare skin about 5 cm below your left nipple.

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Ensure that there is no skin contact between your hand and chest or the recording will not be correct.

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### Method 2: Fingers

Place one finger from your left hand on the lower electrode.

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Place your right index finger on the upper electrodes.

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Ensure that there is no skin contact between your right and left hands or the recording will not be correct.

Press with moderation and steadiness during the recording.

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
The ECG recording can be invalid for the following reasons:

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Do not move / speak during the recording and hold the device still.

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The electrodes must not come into contact with other electrically conducting parts.

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
Do not carry out record your ECG if there is any moisture on the skin like sweat, AFiber showering / bathing.

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Do not use the monitor near strong electromagnetic fields and keep away from radio equipment, mobile phones, telephones and high frequency devices that can impair its functionality.

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If you have a pacemaker or another implanted device, please consult your doctor before recording your own ECG: follow the advice given to you by your doctor

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- Sources:
1. WHO World Health Data Platform „Global Health Estimates: Life expectancy and leading causes of death and disability“, accessed: 2021-01-28, [www.who.int/data/gho/data/themes/mortality-and-global-health-estimates/gh-leading-causes-of-death](http://www.who.int/data/gho/data/themes/mortality-and-global-health-estimates/gh-leading-causes-of-death).
  2. Sunrise Hospital & Medical Center (2017-01-23) „Atrial fibrillation: You may have it and not know it“, accessed: 2021-01-28, [www.lasvegassun.com/native/sunrise/2017/jan/23/atrial-fibrillation-you-may-have-it-and-not-know-it/](http://www.lasvegassun.com/native/sunrise/2017/jan/23/atrial-fibrillation-you-may-have-it-and-not-know-it/).
  3. American Heart Association (2021) „FAQs of Atrial Fibrillation“ accessed: 2021-01-28, [www.heart.org/en/health-topics/atrial-fibrillation/afib-resources-for-patients--professionals/afib-awareness](http://www.heart.org/en/health-topics/atrial-fibrillation/afib-resources-for-patients--professionals/afib-awareness).
  4. Rienstra, Michiel et al. „Symptoms and functional status of patients with atrial fibrillation: state of the art and future research opportunities.“ *Circulation* vol. 125,23 (2012): 2933-43. doi:10.1161/CIRCULATIONAHA.111.069450.
  5. Boriani, Giuseppe et al. „Asymptomatic atrial fibrillation: clinical correlates, management, and outcomes in the EORP-AF Pilot General Registry.“ *The American journal of medicine* vol. 128,5 (2015): 509-18 e2. doi:10.1016/j.amjmed.2014.11.026.
  6. Go, A S et al. „Prevalence of diagnosed atrial fibrillation in adults: national implications for rhythm management and stroke prevention: the AnTicoagulation and Risk Factors in Atrial Fibrillation (ATRIA) Study.“ *JAMA* vol. 285,18 (2001): 2370-5. doi:10.1001/jama.285.18.2370.
  7. Miyasaka, Yoko et al. „Secular trends in incidence of atrial fibrillation in Olmsted County, Minnesota, 1980 to 2000, and implications on the projections for future prevalence.“ *Circulation* vol. 114,2 (2006): 119-25. doi:10.1161/CIRCULATIONAHA.105.595140.
  8. ARNI Institute, T. Balchin (2018) „Atrial Fibrillation raises risk of stroke by 5: Can you tell if you have it?“, accessed: 2021-01-28, [www.arni.uk.com/atrial-fibrillation-raises-risk-of-stroke-by-5-can-you-tell-if-you-have-it/](http://www.arni.uk.com/atrial-fibrillation-raises-risk-of-stroke-by-5-can-you-tell-if-you-have-it/).
  9. American Heart Association (c2019). „AFib Can Happen to Anyone“, accessed: 2021-01-28, [www.heart.org/en/health-topics/atrial-fibrillation/afib-resources-for-patients--professionals/afib-can-happen-to-anyone-infographic](http://www.heart.org/en/health-topics/atrial-fibrillation/afib-resources-for-patients--professionals/afib-can-happen-to-anyone-infographic).

PAUL HARTMANN Middle East  
Dubai Airport Freezone  
P.O. Box 54525 Dubai  
United Arab Emirates  
Phone: +971 4 2996996

me.hartmann.info