
Veroval® compact

BPU 22

EN – BPU 22 Upper Arm Blood Pressure Monitor <i>Instructions for use</i>	2-19
TR – BPU 22 Üst koldan tansiyon ölçme cihazı <i>Kullanım kılavuzu</i>	20-37
BPU 22 جهاز قياس ضغط الدم بالجزء العلوي من الذراع – AE 38-55	تعليمات الاستخدام
HK – BPU 22 上臂式血壓計 <i>使用說明書</i>	56-73
Electromagnetic Compatibility Information	74
Warranty certificate / Garanti belgesi / شهادة الضمان / 保固書	83

HARTMANN

Dear customer

We are pleased that you have chosen to buy a blood pressure monitor from HARTMANN. The Veroval® BPU 22 is a quality product for fully automatic blood pressure measurement on the upper arm of adults and is suitable for clinical and domestic use. Requiring no preliminary settings, this blood pressure monitor easily and automatically inflates for simple, quick and reliable measurement of the systolic and diastolic blood pressure and the pulse rate. The monitor also indicates if it detects an irregular heartbeat.

We wish you all the best for your health.

Please read these instructions for use carefully before using the device the first time because correctly measuring your blood pressure depends on appropriate handling of the device. These instructions for use are designed to teach you each step for self-measurement of blood pressure using the Veroval® BPU 22. You will be given important and useful hints to ensure that you obtain reliable results for your personal blood pressure profile. Use this device in accordance with the information provided in the instructions for use. Please store these instructions carefully and ensure other users can access them. Check that the outside of the package is undamaged and that the contents are complete.

Included in delivery:

- Blood pressure monitor
- Universal upper arm cuff
- 4× 1.5 V AA batteries
- Storage bag
- Instructions for use with warranty certificate


Table of contents

	Page
1. Device and display description	3
2. Important information	4
3. Blood pressure information	9
4. Preparing the measurement	10
5. Measuring your blood pressure	11
6. Memory function	13
7. Explanation of error messages	15
8. Device care	16
9. Warranty conditions	17
10. Contact details for customer queries	17
11. Technical data	18
Warranty certificate	83

1. Device and display description


Blood pressure monitor

- 1 Extra large LCD display
- 2 START/STOP button
- 3 Memory button for User 1
- 4 Memory button for User 2
- 5 Cuff connector socket
- 6 Battery compartment


Cuff

- 7 Cuff with application instructions
- 8 Size indicator to correctly fit cuff
- 9 Cuff plug
- 10 Cuff tube

Display

- 11 Systolic blood pressure
- 12 Diastolic blood pressure
- 13 Pulse rate
- 14 Flashes when the device is measuring and the pulse rate is detected
- 15 Irregular heartbeat
- 16 Display to check the fit of the cuff
- 17 Battery symbol
- 18 Traffic-light system for your values
- 19 Average value (A), morning (AM), evening (PM) / memory allocation number
- 20 User memory
- 21 Date and time display

2. Important information

Signs and symbols

The following symbols are used in the operating instructions, on the packaging and on the type plate of the device and accessories:

Follow instructions for use

Please note

IP21

Protected against solid objects ≥ 12.5 mm and against vertically falling drops of water

Temperature limitation

Air humidity limitation

Protection from electric shock

Dispose of packaging in an environmentally responsible way

Dispose of packaging in an environmentally responsible way

Symbol for the marking of electrical and electronic equipment

Labelling in accordance with Medical Devices Directive 93/42/EEC

Direct current

Manufacturer

Authorised representative in the European Community (EC)

Batch number

Order number

Cardboard recycling code

SN

Serial Number

**Important instructions for use**

- Only use the device to take blood pressure measurements on the upper arm. Do not place the cuff on other areas of the body.
- Use only the supplied cuff or an original replacement cuff. The measured values may otherwise be incorrect.
- Only use the device on persons whose upper arm circumference is within the range indicated for the device.
- If measured values are questionable, repeat the measurement.
- Do not leave the device unattended near small children or persons who cannot operate it themselves. There is a risk of strangulation if the cuff tube becomes wrapped around the neck. Swallowing small parts that have detached from the device can also result in suffocation.
- Do not under any circumstances carry out blood pressure measurements on newborns, babies or toddlers.
- Please do not place the cuff over a wound as this may cause further injuries.
- Do not place the cuff on persons who have had a mastectomy.
- Please note that the pressure build-up in the cuff can lead to temporary disruption of medical devices being used at the same time on the same arm.

- Do not use the blood pressure monitor together with a high frequency surgical device.
- If an intravenous treatment is being carried out or a venous catheter is inserted in the arm, blood pressure measurements can lead to injuries. Never use the cuff on the arm to which these conditions apply.
- The function of the arm on which the cuff is placed may become impaired during inflation.
- If you are carrying out the measurement on another person, please ensure that using the blood pressure monitor does not have a lasting negative effect on the blood circulation.
- Taking measurements too often within a short time period as well as maintaining continuous cuff pressure can disrupt the blood circulation and cause injuries. Please take a break between measurements and do not bend or fold the air tube. If the monitor malfunctions, remove the cuff from the arm.
- Do not use the blood pressure monitor on patients with pre-eclampsia during pregnancy.

Important instructions for self-measurements

- Even minor changes in internal and external factors (e.g. deep breathing, stimulants, talking, excitement, environmental factors) can lead to fluctuations in blood pressure. This explains why deviating values are often measured at the doctor or pharmacist.
- Measurement results depend fundamentally on the measurement location and the position of the patient (sitting, standing, reclining). They are also influenced by physical activity and physiological preconditions of the patient, for example. For comparable values, carry out the measurement in the same location and the same position.
- Diseases of the cardiovascular system can lead to erroneous readings or lower the accuracy of the measurement. The accuracy may also be affected if you have very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or tremors.

Consult your doctor before measuring your own blood pressure if ...

- you are pregnant. Blood pressure may vary during pregnancy. Regular blood pressure monitoring is particularly important if you have high blood pressure because elevated values may affect foetal development. Check with your doctor in any case to find out if and when you should measure your own blood pressure, especially if you suffer from pre-eclampsia.

- you have diabetes, liver function disorders or narrowing of blood vessels (e.g. arteriosclerosis, peripheral arterial occlusive disease) because deviating measurements may occur in such cases.
- you have certain blood disorders (e.g. haemophilia) or severe circulatory disorders or are taking blood-thinning drugs.
- you wear a cardiac pacemaker because it may lead to deviating measured values. The blood pressure monitor itself has no impact on a cardiac pacemaker. Please note that the pulse rate displayed is not suitable for checking the frequency of cardiac pacemakers.
- you tend to bruise easily and/or are sensitive to pressure pain.
- you suffer from severe irregular heartbeat or cardiac arrhythmias. Because of the oscillometric measurement method that is used, in some cases incorrect measured values are obtained or the blood pressure cannot be measured.
- If this symbol  appears often, it may indicate heart rhythm disorders (arrhythmias). If this is the case, please consult your doctor. Severe heart rhythm disorders may sometimes produce false measurements or lower the measurement accuracy. Please discuss with your doctor whether taking your own blood pressure readings is suitable for you.
- Your self-measured values are for your information only – they are not a substitute for a medical examination! Discuss your measured values with your doctor and under no circumstances make your own medical decisions based on these measurements (e.g. medications and their dosages)!

- Measuring your own blood pressure is not a substitute for medical treatment! Do not interpret your measured values by yourself and do not use them for self-prescribed treatment. Take measurements based on the instructions of your doctor and trust his or her diagnosis. Take medications as prescribed by your doctor and never change the dose on your own. Discuss the appropriate time to measure your blood pressure with your doctor.

Handling batteries

- Observe the polarity labels plus (+) and minus (-).
- Only use high-quality batteries (see specification in Chapter 11 'Technical data'). If you use low-quality batteries, we cannot guarantee the specified number of measurements.
- Never mix old and new batteries or batteries from different manufacturers.
- Remove empty batteries immediately.
- Replace batteries if the battery symbol  remains illuminated.
- Always replace all the batteries at the same time.
- If the device remains unused for a longer period of time, you should remove the batteries to prevent possible leakage.

Batteries

■ Choking hazard

Small children could swallow batteries and suffocate on them. Keep batteries out of the reach of children!

■ Risk of explosion

Do not throw batteries into a fire.

■ Batteries must not be charged or short-circuited.

■ If a battery has leaked, wear protective gloves and clean the battery compartment with a dry cloth. If liquid from a battery cell comes into contact with skin or eyes, clean the affected area with water and seek medical attention if necessary.

■ Protect batteries from excessive heat.

■ Do not disassemble, open or crush batteries.

Safety information concerning the device

- This blood pressure device is not waterproof!
- This blood pressure device is made of high-quality electronic precision components. The accuracy of the measured values and the lifetime of the device depend on careful handling.
- Protect the device from strong shocks, impacts or vibrations and do not let it drop on the floor.
- Do not excessively bend or kink the cuff or the air tube.

- Never open the device. The device must not be modified, dismantled or repaired by the user. Repairs may only be carried out by authorised specialists.
- Never inflate the cuff when it is not properly applied to the upper arm.
- Use the device only with the approved upper arm cuff. The device may otherwise sustain external or internal damage.
- The cuff tube may only be removed from the unit by pulling on the corresponding connector plug. Never pull on the tube itself!
- Do not expose the device to extreme temperatures, humidity, dust or direct sunlight because this may cause it to malfunction.
- Keep the packaging, batteries and device out of the reach of children.
- Please comply with the storage and operating conditions defined in Chapter 11 'Technical data'. Storage or use outside the specified temperature and humidity range can affect the accuracy of the measurement or the function of the device.
- Do not use the device near strong electromagnetic fields and keep away from radio equipment or mobile phones. Portable and mobile high-frequency and communication devices, such as telephones and mobile phones, can impair the functionality of this electronic medical device.

Instructions for the measurement function check

Every Veroval® device has been carefully tested by HARTMANN for measurement accuracy and has been developed with a view to a long

service life. We recommend carrying out a metrological check every 2 years for **devices in professional use**, for example, in pharmacies, medical practices or hospitals. Please also observe the national regulations specified by the legislator. Metrological checks may only be carried out by the competent authorities or authorised maintenance providers against cost reimbursement.

Disposal information

- To protect our environment, empty batteries must not be disposed of in household waste. Please comply with the relevant waste disposal regulations or use public collection points.
- This product is subject to the European Directive 2012/19/EU on Used Electrical and Electronic Equipment and is marked accordingly. Never dispose of electronic equipment in your household waste. Please learn about your local regulations on the proper disposal of electrical and electronic products. Proper disposal protects the environment and human health.

3. Blood pressure information

To determine your blood pressure you need to measure two values:

- **SYS** – Systolic (upper) blood pressure, which is produced when the heart contracts and pumps blood into the blood vessels.
- **DIA** – Diastolic (lower) blood pressure, which is the value measured when the heart muscle is dilated and fills again with blood.
- Blood pressure values are displayed in mmHg.

So that you can understand your results more easily, a traffic-light system is available on the left side of the Veroyal® BPU 22 that directly indicates the result, making it easier to categorise the measured value. The World Health Organisation (WHO) and the International Society of Hypertension (ISH) have developed the following summary for classifying blood pressure values:

Results indicator	Evaluation	Systolic pressure	Diastolic pressure	Recommendation
red	Grade 3 hypertension	above 179 mmHg	above 109 mmHg	Consult a doctor
orange	Grade 2 hypertension	160 - 179 mmHg	100 - 109 mmHg	
yellow	Grade 1 hypertension	140 - 159 mmHg	90 - 99 mmHg	Regular examinations by a doctor
green	Normal limit values	130 - 139 mmHg	85 - 89 mmHg	
green	Normal	120 - 129 mmHg	80 - 84 mmHg	Self-assessment
green	optimal	up to 119 mmHg	up to 79 mmHg	

Source: WHO, 1999 (World Health Organisation)

4. Preparing the measurement

Inserting / changing the batteries

■ Open the battery cover on the underside of the device. Insert the batteries (see Chapter 11 'Technical data'). Ensure correct polarity ('+' and '-') when inserting batteries. Close the battery lid. 12 h or 24 h flashes on the display. Set date and time as described below.

■ If the 'Change battery' symbol  is permanently illuminated, blood pressure can no longer be measured and you must replace all the batteries.

Setting the date and time



Be sure to set the date and time correctly. This is the only way to save your measured values correctly with date and time for subsequent retrieval.

■ To switch to setting mode, reinsert the batteries or hold the START/STOP button  down for 5 seconds. Proceed as follows:

Hours



The hour format flashes on the display.

- Select the desired hour format using the Save buttons  /  and confirm using the START / STOP button .

Date

The display shows the year (a), the month (b) and the day (c).



- Depending on the display, use the Save buttons  /  to select the year, the month or the day and confirm using the START / STOP button .

If the 12-hour format is selected, the month is displayed before the day.

Time

The hour (d) and the minutes (e) flash consecutively on the display.

- Depending on the display, use the Save buttons  /  to select the current hour or minutes and confirm using the START / STOP button .

The device switches off automatically once all data has been set.

5. Measuring your blood pressure

Application of the cuff

- Before applying the cuff, insert the cuff socket into the connector socket on the left side of the device.
- Do not mechanically constrict, compress or bend the sleeve tube.
- The blood pressure must be measured on a bare upper arm. If the cuff is completely open, thread the end of the cuff through the metal ring to form a loop (see Fig. 1). Position the cuff on the upper arm so that the lower edge is 2-3 cm above the crook of the arm and above the artery (see Fig. 2). The tube points to the centre of the palm.


Fig. 1


Fig. 2

- Now bend your arm slightly, grasp the free end of the cuff, wrap it firmly round your arm and close the hook-and-loop fastener.

The cuff should fit securely but not tightly. You should be able to insert two fingers between the arm and the cuff. Make sure the tube is not bent or damaged.

Important: It is essential that the cuff is correctly applied to obtain correct measurements. The cuff included with the device is a universal cuff for upper arm circumferences of 22 cm to 42 cm. The white arrow must point to an area within the sizing scale. If it is outside the range, it is not possible to guarantee a correct measurement result.

This innovative Veroval® device with Comfort Air Technology ensures comfortable measurement. Inflation to 190 mmHg is required for the first measurement. For subsequent measurements, inflation pressure is adjusted individually based on previously measured blood pressure values. This allows for more comfortable upper arm measurements.

Carrying out the measurement

- Rest for approx. 5 minutes before measurement.
- You should take your blood pressure in a quiet place, in a relaxed and comfortable seated position.
- Measurement can be taken on the right or left arm. We recommend you perform the measurement on the left upper arm. Over the long term, the arm giving higher readings should be used for blood pressure monitoring. However, if there is a very clear

difference between readings on either arm, you need to check with your doctor which arm you should use for the measurement.

- Always measure on the same arm and place your forearm in a relaxed position on a support.
- We recommend that you measure your blood pressure while sitting with your back supported by the back of the chair. Place both feet flat on the floor next to each other. Do not cross your legs. Relax your forearm and hand with your palm facing upwards on a support and make sure the cuff is at the height of your heart.
- Always take measurements at the same time of day. Only regular measurements taken at the same time each day over a longer period of time give a meaningful assessment of blood pressure values.
- Do not measure when you have a strong urge to urinate. A full bladder can increase blood pressure by about 10 mmHg.
- Do not take your blood pressure after taking a bath or playing sport.
- Do not eat, drink or exercise for at least 30 minutes before the measurement.
- Please wait for at least one minute between two measurements.
- Start a measurement only after applying the cuff. Press the START/STOP button . The appearance of all display segments followed by the time and date indicates that the device is carrying out its self-test and is ready for use.
- Check the display segments for completion (see Chapter 1).
- After about 0,5 seconds the cuff will automatically inflate. If this inflation pressure is insufficient or if the measurement is inter-


rupted, the device continues to pump in increments of 40 mmHg until a high enough pressure is reached. During inflation, the results indicator on the left also increases at the same time.

- If the cuff has been applied securely enough on the arm, the cuff symbol appears  on the display. If the cuff symbol does not appear in the display, the cuff has not been fitted securely enough and the error message 'E3' appears on the display after a few seconds.

Important: Please do not speak or move during the measurement.

- As the air pressure in the cuff decreases, the heart symbol  and the dropping pressure in the cuff are shown on the display.
- After the measurement is complete, the systolic and diastolic blood pressure as well as the pulse rate (see Fig.) appear simultaneously.

- In addition to the measured values, the time, date, the associated user memory [1] or [2] as well as the associated memory numbers (e.g. (M) 05) appear. The measured value is automatically stored in the displayed user memory. As long as the result of the measurement is displayed, you can assign values to the corresponding user memory by pressing the [1] or [2] button. If you do not assign them, the measured value is automatically stored in the displayed user memory. You can use the results indicator on the left of the display to classify your measurement result (see table in Chapter 3 'Blood pressure information').
- To switch off the device, press the START/STOP button (1), otherwise the device will automatically switch off after 1 minute.

If you wish to stop the measurement for any reason, simply press the START/STOP button. (1) The cuff inflation or measurement is interrupted and the pressure is automatically released.

- If this symbol  is shown in the display, the unit has detected an irregular heartbeat during the measurement. However, the measurement may also have been disrupted by body movement or speaking. It is best to repeat the measurement. If you see this symbol regularly when measuring your blood pressure, we recommend that your doctor checks your heart rhythm.

6. Memory function

User memory

- The Veroval® BPU 22 stores up to 100 measurements in each user memory. Once all memory slots are occupied, the oldest value is deleted.
- Memory recall is activated by pressing the [1] or [2] button when the device is switched off. For values saved in the first user memory, press the [1] button, for the second user memory press the [2] button.

Average values

- After selecting the particular user memory, the corresponding symbol [1] or [2] and an \bar{M} appear on the display. The average value of all the data stored for the particular user memory is displayed (see Fig. 1).

Fig. 1

- By pressing the **[1]** button again (or **[2]** button if you are in user memory 2), the average values for all morning measurements 'AM' (5 am to 9 am) taken on the last 7 days appear.

- By pressing the **[1]** button again (or **[2]** button if you are in user memory 2) the average values for all evening measurements 'PM' (6.00 pm to 8.00 pm) taken on the last 7 days appear.

Individual measured values

- By pressing the **[1]** button again (or **[2]** button, if you are in user memory 2), all memory values can be retrieved one after the other, starting with the most recent measured value.

- If an irregular heartbeat was detected during measurement, that information  is also stored and displayed when the measured value is retrieved from the device's memory along with the systolic and diastolic blood pressure readings, pulse rate, time, date and year.
- You can cancel the memory recall at any time by pressing the START/STOP button **[1]**. Otherwise, the device will switch off automatically after a few seconds.
- Even if the power supply is interrupted, e.g. when changing the batteries, the stored values will still be available.

Deleting saved readings

All data saved for each user can be separately deleted for user memory **[1]** and user memory **[2]**. To do so, press the button of the corresponding user memory (**[1]** or **[2]**). The average value will then appear on the display.

Press and hold the user memory button for 5 seconds. 'CL 00' appears on the display. All data in the selected user memory are now deleted. If you release the button ahead of time, no data will be deleted.

7. Explanation of error messages

Error that has occurred	Possible causes	Remedy
Monitor cannot be switched on	Batteries are missing, inserted incorrectly or flat.	Check batteries and insert four identical new batteries if necessary.
Cuff will not inflate	Cuff connector plug is incorrectly positioned in the connector socket of the monitor.	Check the connection between the cuff connector plug and socket.
	Wrong cuff type connected.	Verify that only an approved Veroval® cuff and connector are being used.
E1	Pulse could not be detected correctly.	Check whether the cuff has been applied correctly. Do not talk or move during the measuring procedure.
E2	The pulse could not be detected.	Check whether the cuff has been positioned correctly. Do not talk or move during the measuring procedure.
E3	The cuff was applied too firmly or too loosely.	Apply the cuff so that two fingers fit between the cuff and the upper arm.
		Air tube is not correctly connected to the monitor. Check that the connector plug is correctly positioned. If this error occurs frequently, you should use a new cuff.

Error that has occurred	Possible causes	Remedy
E4 / E6	There is a system error.	Contact customer service if this error message appears.
E5	Inflation pressure is higher than 300 mmHg.	Please rest for at least 1 minute and take the measurement again.
	The batteries are almost empty.	Replace batteries.
Implausible measured values	Implausible measured values often occur due to inappropriate handling of the device or mistakes during the measuring procedure.	Please observe Chapter 5 'Measuring your blood pressure' as well as the safety information. Then repeat the measurement.

Switch the device off if an error symbol appears. Check for all possible causes and note the instructions for self-measurement from Chapter 2 'Important information'. Relax for a minute and then take the measurement again.

8. Device care

- Only ever use a soft, damp cloth to clean the device. Please do not use thinner, alcohol, detergents or solvents.
- The cuff can be cleaned carefully using a slightly damp cloth and mild soap solution. Do not completely immerse the cuff in water.

- It is recommended to clean and disinfect the cuff regularly or after each use, especially when used by several users, to prevent infection. The cuff should be disinfected, particularly on the inside, by wiping with a disinfectant. Use a disinfectant that is compatible with the cuff materials, e.g. 75% ethanol or isopropyl alcohol.

To protect from external influences, keep the device, the cuff and these instructions in the storage bag.

- If you store a device and cuff, do not place any heavy objects on the device and the cuff. Remove the batteries.

9. Warranty conditions

- We offer a 3-year warranty for this premium quality blood pressure monitor from the day of purchase and in accordance with the conditions listed below.
- Claims must be made during the warranty period. The date of purchase must be documented by an appropriately completed and stamped warranty certificate or proof of purchase.
- Within the warranty period, HARTMANN will replace or repair free of charge any faulty device components caused by material faults or manufacturing errors. This does not extend the warranty period.
- The device is intended for the purpose described in this instruction manual only.
- The warranty does not cover damage resulting from improper use or unauthorised interference. Accessory parts that are subject to wear and tear (batteries, cuffs, etc.) are excluded from the warranty. Claims for compensation are limited to the value of the merchandise; compensation for consequential damage or injury is expressly excluded.

- In warranty cases, please send the device with cuff together with the fully completed and stamped warranty certificate or purchase receipt directly or via your dealer to the Customer Services department for your country.

10. Contact details for customer queries

ZA - HARTMANN South Africa
Epsom Avenue, Northriding,
2169 Johannesburg
www.veroval.info
Tel. +27 860 4278 6266
phzahelpdesk@hartmann.info

11. Technical data

Model:	Veroval® BPU 22
Type:	GCE606
Measuring method:	oscillometric
Display range:	0 - 300 mmHg
Measuring range:	Systolic (SYS): 50–280 mmHg Diastolic (DIA): 30–200 mmHg Pulse: 40–199 pulse beats per minute Displaying of correct values outside the measuring range cannot be guaranteed.
Display unit:	1 mmHg
Technical measuring accuracy:	Cuff pressure: ± 3 mmHg Pulse: $\pm 5\%$ of the displayed pulse rate
Clinical measuring accuracy:	complies with the requirements of DIN EN 1060-4; Korotkoff validation method: Phase I (SYS), Phase V (DIA)
Operating mode:	continuous mode
Nominal voltage:	DC 6V

Power supply:	4× 1.5 V alkaline manganese (AA/LR06) batteries
Expected service life:	20,000 measurements
Battery capacity:	about 1,000 measurements
Protection against electric shock:	Medical-electrical equipment with internal power supply (only when batteries are used); applied part: type BF
Protection against harmful ingress of water or solid materials:	IP21 (not protected against moisture)
Inflation pressure:	about 190 mmHg for first measurement
Automatic switch-off function:	1 minute after measurement / otherwise 30 sec.
Cuff:	Veroval® cuff for BPU 22, cuff for arm circumference of 22 to 42 cm
Memory capacity:	2 × 100 measurements with average of all measurements and morning / evening average of the last 7 days

Operating conditions:	Ambient temperature: +10 °C to +40 °C Relative humidity: < 90 %, non-condensing Air pressure: 800 hPa – 1050 hPa
Storage / transport conditions:	Ambient temperature: –20°C to +55°C Relative humidity: <90 %, non-condensing
Serial number:	inside the battery compartment
Reference to standards:	IEC 60601-1; IEC 60601-1-2
Weight:	approx. 220 g (without batteries)
Dimensions:	approx. 134(L) x 48(W) x 91(H) mm

Legal requirements and guidelines

- The Veroval® BPU 22 conforms to the European regulations, which are based on the European Medical Device Directive 93/42/EEC, and bears the CE mark.
- The device complies, for example, with the European standard EN 1060: Non-invasive sphygmomanometers – Part 3: Additional requirements for electro-mechanical blood pressure measuring systems and the standard IEC 80601-2-30.
- Clinical testing of measuring accuracy was performed in accordance with the European Standard EN 1060-4.
- Beyond statutory requirements, the device has been clinically validated by the ESH-IP2 Protocol of the European Society of Hypertension (ESH).

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Electromagnetic compatibility

Table 1

For all medical-electrical (ME) DEVICES and ME SYSTEMS

Guidance and manufacturer's declaration - electromagnetic emissions

The Veroval® BPU 22 is intended for use in the electromagnetic environment specified below. The customer or the user of the Veroval® BPU 22 should assure that it is used in such an environment

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Veroval® BPU 22 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class [B]	The Veroval® BPU 22 is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Table 2**Guidance and manufacturer's declaration - electromagnetic Immunity**

The Veroval® BPU 22 is intended for use in the electromagnetic environment specified below. The customer or the user of the Veroval® BPU 22 should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV input/output lines: ±1 kV	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV. line(s) to earth: ±2 kV. 100 kHz repetition frequency	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycle 70% U_T (30% dip in U_T) for 25 cycle <5% U_T (>95% dip in U_T) for 5s	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Table 3**Guidance and manufacturer's declaration - electromagnetic Immunity**

The Veroval® BPU 22 is intended for use in the electromagnetic environment specified below. The customer or the user of the Veroval® BPU 22 should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance	
Conducted RF IEC61000-4-6	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the Veroval® BPU 22, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distances: $d = 1.2 \sqrt{P}$; $d = 2 \sqrt{P}$	
Radiated RF IEC61000-4-3	10V/m	10V/m	80 MHz to 800 MHz: $d = 1.2 \sqrt{P}$ 800 MHz to 2.7 GHz: $d = 2.3 \sqrt{P}$	Where, P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, d is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.

To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Veroyal[®] BPU 22 is used exceeds the applicable RF compliance level above, the Veroyal[®] BPU 22 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Veroyal[®] BPU 22.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4**Recommended separation distances between portable and mobile RF communications equipment and the Veroval® BPU 22**

The Veroval® BPU 22 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Veroval® BPU 22 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Veroval® BPU 22 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz (out ISM and amateur radio bands) $d = 1.2 \sqrt{P}$	150 kHz to 80 MHz (in ISM and amateur radio bands) $d = 2 \sqrt{P}$	80MHz to 800MHz $d = 1.2 \sqrt{P}$	800MHz to 2.7GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.2	0.12	0.23
0.1	0.38	0.632	0.38	0.73
1	1.2	2	1.2	2.3
10	3.8	6.32	3.8	7.3
100	12	20	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Table 5**Guidance and manufacturer's declaration - electromagnetic Immunity**

The Veroval® BPU 22 is intended for use in the electromagnetic environment specified below. The customer or the user of the Veroval® BPU 22 should assure that it is used in such an environment.

	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 –390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27
	450	430 –470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation b) 217 Hz	0,2	0,3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0,3	28
	870						
	930						
	1720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0,3	28
	1845						
	1970						

	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	2450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0,3	28
	5240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0,2	0,3	9
	5240						
	5785						

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

- a) For some services, only the uplink frequencies are included.
- b) The carrier shall be modulated using a 50 % duty cycle square wave signal.
- c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E = \left[\frac{6}{d} \right] \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

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Dealer's stamp · Yetkili satıcı kaşesi · ختم التاجر · 專賣店印章

Date of revision: 2017-10

Made under the control of PAUL HARTMANN AG

Globalcare Medical Technology Co., Ltd.
7th Building, 39 Middle Industrial Main Road
European Industrial Zone, Xiaolan Town
528415 Zhongshan City, Guangdong Province
P.R.C.

Donawa Lifescience Consulting Srl
Piazza Albania, 10
00153 Rome
Italy

758.989-0618