# Veroval®

ECG and blood pressure monitor

EN – ECG and blood pressure monitor

*Instructions for use*..... 2 - 34

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



HARTMANN

# 1. Description of the device and display

## **ECG and blood pressure monitor**

- 1 Extra large LCD display
- 2 SETTINGS button
- 3 START/STOP button
- 4 SAVE button
- 5 Cuff connector socket
- 6 Upper electrodes for the ECG measurement
- 7 Lower electrode with detection button for ECG measurement
- 8 Battery compartment
- 9 USB *Bluetooth*® adapter

## **Cuff**

- 10 Secure fit cuff (a) with application instructions (b)
- 11 Grip tab to pull on the cuff
- 12 Size indicator to correctly fit cuff
- 13 Cuff plug
- 14 Cuff tube
- 15 Lower arm opening for correct and secure fit

## **Display**

- 16 Traffic-light system for your blood pressure readings
- 17 *Bluetooth*® display, indicates the *Bluetooth*® connection status
- 18 User memory
- 19 Data transfer symbol, flashes when data is being transferred
- 20 Alarm indicator, shows whether the alarm is activated or deactivated
- 21 During the ECG measurement: 30-second countdown.  
While the results are being displayed: memory position number and month/  
day/hour/minute.
- 22 Battery symbol
- 23 Measurement mode: ECG (ECG record) or BP (blood pressure measurement)
- 24 ECG analysis
- 25 Pulse rate
- 26 Heart symbol, flashes during the ECG recording when the device is measuring  
and the pulse is detected
- 27 Diastolic blood pressure
- 28 Systolic blood pressure

## 2. Important notes

### Signs and symbols

	Follow instructions for use
	Please note
<b>IP22</b>	Protected against penetration by solid objects with a diameter $\geq 12.5$ mm. Protected against water drops falling at an angle of up to $15^\circ$ to the monitor.
	Temperature limitation
	Humidity limitation
	Protection from electric shock
	Dispose of in an environmentally responsible way
	Dispose of in an environmentally responsible way
	Symbol regarding electrical and electronic equipment labelling
<b>CE</b> <sub>0120</sub>	Labelling according to Directive 93/42/EEC on Medical Devices
	Direct current
	Manufacturer
<b>EC REP</b>	Authorised representative in the European Community (EC)
<b>LOT</b>	Batch number
<b>REF</b>	Order number
	Cardboard recycling code

**SN** Serial number

## **Important information about use**

- Only use the device to take blood pressure measurements on the upper arm or for recording an ECG using the methods described in Chapter 7: Recording an ECG. Do not place the cuff on other areas of the body.
- Use only the supplied cuff or an original replacement cuff. The measured values may otherwise be incorrect.
- Only use the device to measure the blood pressure of persons whose upper arm circumference is within the range indicated for the device.
- At an ambient temperature of 40°C both the cuff and the ECG electrodes reach a temperature of 42°C.
- If measured values are questionable, repeat the measurement.

- Do not leave the monitor unattended near small children or persons who cannot operate it themselves. There is a risk of strangulation if the cuff tube becomes wrapped around the neck.  
Swallowing small parts can result in suffocation.
- Do not under any circumstances carry out ECG recording or blood pressure measurements on newborns, babies or toddlers.
- Please do not place the cuff over a wound as this may cause further injuries.
- Do not place the cuff on people who have had a mastectomy.
- Please note that the pressure build-up in the cuff can lead to temporary disruption to medical devices being used at the same time on the same arm.
- Do not use the ECG and blood pressure monitor together with a high frequency surgical device.
- If an intravenous treatment is being carried out or a venous catheter is inserted in the arm, blood pressure measurements can lead to injury. Never use the cuff on the arm to which these conditions apply.

- The function of the arm on which the cuff is placed may become impaired during inflation.
- If you are carrying out the measurement on another person, please ensure that using the ECG and blood pressure monitor does not have a lasting negative effect on the blood circulation.
- Taking measurements too often within a short time period as well as maintaining continuous cuff pressure can disrupt the blood circulation and cause injury. Please take a break between measurements and do not bend or fold the air tube. If the monitor malfunctions, remove the cuff from the arm.
- Do not use the ECG and blood pressure monitor on patients with pre-eclampsia during pregnancy.

### **Important instructions for measuring your own blood pressure**

- Even minor changes in internal and external factors (e.g. deep breathing, stimulants, talking, excitement, environmental factors) can lead to fluctuations in blood pressure. This explains why anomalous values are often measured by a doctor or pharmacist.
- Measurement results depend fundamentally on where measurement location and the position of the patient (sitting, standing, reclining). They are also influenced by physical activity and physiological preconditions of the patient, for example. For comparable values, carry out the measurement in the same location and the same position.
- Diseases of the cardiovascular system can lead to erroneous readings or lower the accuracy of the measurement. The accuracy may also be affected if you have very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or tremors.

### **Consult your doctor before measuring your own blood pressure if ...**

- you are pregnant. Blood pressure may vary during pregnancy. Regular blood pressure monitoring is particularly important if you have high blood pressure because elevated values may affect foetal development. Check with your doctor in any case to find out if and when you should measure your own blood pressure, especially if you suffer from pre-eclampsia.

- you suffer from diabetes, liver disorders or congestion of blood vessels (e.g. arteriosclerosis, peripheral arterial occlusive disease) because anomalous measured values may occur in such cases.
- you suffer from certain blood diseases (e.g. haemophilia) or severely impaired blood flow or if you take blood-thinning drugs.
- you wear a cardiac pacemaker because it may lead to anomalous measured values. Please note that the pulse rate displayed is not suitable for checking the frequency of cardiac pacemakers.
- you tend to bruise easily and/or are sensitive to pressure pain.
- you suffer from severe irregular heartbeat or cardiac arrhythmias. Because of the oscillometric measurement method that is used, in some cases incorrect measured values are obtained or the blood pressure cannot be measured.
- The values measured by you are for your information only – they are not a substitute for a medical examination! Discuss your measured values with your doctor and under no circumstances should you make personal medical decisions based on these measurements (e.g. medications and their dosages)!
- Measuring your own blood pressure is not a substitute for medical treatment! Do not interpret your measured values by yourself and do not use them for self-prescribed treatment. Take measurements based on the instructions of your doctor and trust his or her diagnosis. Take medications as prescribed by your doctor and never change the dose on your own. Discuss the appropriate time to measure your blood pressure with your doctor.

### **Important instructions for recording your own ECG**

- The device displays changes in the heart rhythm. These may have a number of different causes which may be harmless but which may also be triggered by diseases of varying severity. If you suspect you have a medical condition, please seek medical advice.
- Do not carry out any measurements if there is any moisture on the skin such as sweat or after showering/bathing.
- If you have a pacemaker or another implanted device, please consult your doctor before recording your own ECG: follow the advice given to you by your doctor.
- The recorded electrocardiograms reflect the function of the heart at the time of the measurement. Previous or subsequent changes therefore cannot necessarily be detected.

- The ECG measurements carried out using the device cannot identify every single heart disease. Immediately seek medical treatment, regardless of the measurement result, if you notice any symptoms that may indicate an acute heart disease. If you should notice any of the following symptoms, please have them clarified medically (this list is not exhaustive): pain on the left side or a feeling of pressure in the chest or abdominal region; radiating pain in the mouth/jaw/face, shoulder, arm or hand; pain around the back; nausea; burning in the chest area; tendency to collapse; difficulty breathing; racing heart or irregular heart rhythm.
- Do not use any defibrillators while recording an ECG.
- Do not record an ECG while undergoing an MRI examination.
- Recording your own ECG is not a substitute for medical treatment! Do not interpret the electrocardiograms by yourself and do not use them for self-prescribed treatment. Take measurements based on the instructions of your doctor and trust his or her diagnosis. Take medications as prescribed by your doctor and never change the dose on your own.
- The device is not a substitute for either a medical examination of your heart function, nor recording of a medical electrocardiogram that must be obtained using more elaborate measuring equipment.
- The device does not diagnose any medical conditions that may be the possible cause of the change observed in the ECG. This can only be done by your treating doctor.
- It is recommended that you record the ECG curves measured and give them to your treating doctor if necessary.

### **Power supply (batteries)**

- Observe the polarity labels plus (+) and minus (-).
- Only use high-quality batteries (see specification in Chapter 14 Technical data). If you use low-quality batteries, we can no longer guarantee the specified number of measurements.
- Never mix old and new batteries or batteries from different manufacturers.
- Remove empty batteries immediately.
- Replace batteries if the battery symbol  remains illuminated.
- Always replace all the batteries at the same time.
- If you do not intend to use the monitor for a longer period of time, you should remove the batteries to prevent possible leakage.

## Batteries

### ■ Choking hazard

Small children could swallow batteries and suffocate on them. Keep the packaging, batteries and the monitor out of the reach of children.

### ■ Risk of explosion

Do not throw batteries into a fire.

■ Batteries must not be charged or short-circuited.

■ If a battery has leaked, wear protective gloves and clean the battery compartment with a dry cloth. If liquid from a battery cell comes into contact with skin or eyes, clean the affected area with water and seek medical attention if necessary.

■ Protect batteries from excessive heat.

■ Do not disassemble, open or crush batteries.

## Safety information concerning the monitor

■ This ECG and blood pressure monitor is not waterproof!

■ This ECG and blood pressure monitor is made of high-quality electronic precision components. The accuracy of the measured values and the lifetime of the monitor depend on careful handling.

■ Protect the monitor from violent impact, knocks or vibrations and do not let it drop on the floor.

■ Do not excessively bend or fold the cuff or the air tube.

■ Never open the monitor. You may neither modify nor repair the device yourself. Repairs may only be carried out by authorised specialists.

■ Never inflate the cuff when it is not properly applied to the upper arm.

■ Use the monitor only with the approved upper arm cuff. The monitor may otherwise sustain external or internal damage.

■ The cuff tube may only be removed from the unit by pulling on the corresponding connector plug. Never pull on the tube itself!

■ Do not use the device on persons with sensitive, irritated or injured skin or allergies.

■ This device is not approved for using in children who weigh less than 10 kg.

- The electrodes on the device must not come into contact with other electrically conducting parts (including earth wires).
- The device must not be exposed to static electricity. Always ensure that you have not built up any static electricity before operating the device.
- Do not expose the monitor to extreme temperatures, humidity, dust or direct sunlight because this may cause it to malfunction.
- Keep the packaging, batteries and monitor out of the reach of children.
- Please comply with the storage and operating conditions defined in Chapter 14 Technical data. Storage or use outside the specified temperature and humidity range can affect the accuracy of the measurement or the function of the monitor.
- Do not use the monitor near strong electromagnetic fields and keep away from radio equipment or mobile phones. Portable and mobile high-frequency and communication devices, such as telephones and mobile phones, can impair the functionality of this electronic medical device.

### Instructions for the metrological check

Every Veroval® device has been carefully tested by HARTMANN for measurement accuracy and has been developed with a view to a long service life. We recommend carrying out a metrological check every 2 years for **devices in professional use**, for example, in pharmacies, medical practices or hospitals. Please also observe the national regulations specified by the legislator. Metrological checks may only be carried out by the competent authorities or authorised maintenance providers against cost reimbursement.

### Instructions for calibration mode for the blood pressure function

A functional check of the monitor can be carried out on people or using a suitable simulator. A metrological check involves testing for leaks in the pressure system and possible deviations in the pressure reading. To reach the calibration mode, hold the START/STOP button  for 6 seconds. At the same time, press the Settings button. Release the button after a few seconds and wait 10 seconds. '

10

### Disposal information

- To protect our environment, empty batteries must not be disposed of in household waste. Please comply with the relevant waste disposal regulations or use public collection points.
- This product is subject to the European Directive 2012/19/EU on Used Electrical and Electronic Equipment and is marked accordingly. Never dispose of electronic equipment in your household waste. Please obtain information about local regulations on the proper disposal of electrical and electronic products. Proper disposal protects the environment and human health.

## 3. Blood pressure information

To determine your blood pressure you need to measure two values:

- The systolic (upper) blood pressure: this is produced when the heart contracts and pumps blood into the blood vessels.
- The diastolic (lower) blood pressure: this is the value measured when the heart muscle is dilated and fills again with blood.
- Blood pressure values are displayed in mmHg.

In order to easily assess your results, an indicator can be seen on the right side of the display of the Veroval® ECG and blood pressure monitor that directly shows the result, making it easier to categorise the measured value. The World Health Organization (WHO) and the International Blood Pressure Society (ISH) have developed the following summary for the classification of blood pressure values:

Results indicator	Evaluation	Systolic pressure	Diastolic pressure	Recommendation
<b>3</b>	Grade 3 hypertension	above 180 mmHg	above 110 mmHg	Consult a doctor
<b>2</b>	Grade 2 hypertension	160–179 mmHg	100–109 mmHg	
<b>1</b>	Grade 1 hypertension	140–159 mmHg	90–99 mmHg	Regular examinations by a doctor
<b>N</b>	Normal limit values	130–139 mmHg	85–89 mmHg	
<b>N</b>	Normal	120–129 mmHg	80–84 mmHg	Self-assessment
<b>N</b>	Optimal	up to 119 mmHg	up to 79 mmHg	

Source: WHO, 1999 (World Health Organization)

- Established hypertension (high blood pressure) is defined as having a systolic value of at least 140 mmHg and/or a diastolic value of at least 90 mmHg.
- In general, low blood pressure (hypotension) is defined as a systolic value below 105 mmHg and a diastolic value of 60 mmHg. However, the threshold between normal and low blood pressure (hypotension) is not as accurately specified as the threshold for high blood pressure (hypertension). Hypotension may be associated with symptoms such as dizziness, fatigue, tendency to faint, visual disturbances or high pulse rate. In order to ensure that hypotension or the corresponding symptoms are not attendant symptoms of serious illness, a doctor should be consulted in case of doubt.

Constantly elevated blood pressure significantly increases the risk of developing other medical conditions. The physical consequences of high blood pressure, e.g. heart attack, stroke and organ damage, are the most common causes of death worldwide. Daily blood pressure monitoring is therefore an important measure that can help to protect you from these risks. You should discuss this with your doctor, especially if you have frequently elevated or borderline blood pressure values (see page 49). (The Veroval® medi.connect software lets you share your readings easily with your doctor by email or as a print-out – see Chapter 9 Transferring measured values to Veroval® medi.connect). Your doctor will initiate appropriate measures.

## 4. Information about the ECG

Excitatory processes in the heart can be assessed using an electrocardiogram. A contraction of the heart muscle is always caused by an electrical stimulus. The resulting changes in the voltage can be detected on the surface of the body and their course can be recorded and graphed.

The Veroval® ECG and blood pressure monitor is a single-channel ECG that can detect different arrhythmias.

The basic recording and assessment is carried out by the Veroval® ECG and blood pressure monitor. The device automatically evaluates whether the heart rhythm is too fast (Fig. 3), too slow (Fig. 2), irregular (Fig. 5 and 6), associated with pauses (Fig. 4) or normal (Fig. 1). You will see a clear symbol on the display showing the evaluation immediately after the measurement (more information about the symbols used can be found in Chapter 7 Recording an ECG).

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

Fig. 6

Using the Veroval® medi.connect software, the results recorded can be graphed and printed out for you to present to your doctor. This documentation of patients' results is suitable for medical practitioners and allows your doctor to make a rapid further diagnosis.

## 5. Preparing to measure your blood pressure and record an ECG

Remove the protective film on the display and the electrodes before the first measurement.

### Inserting / changing the batteries

- Open the battery cover on the underside of the device (see Fig. 1). Insert the batteries (see Chapter 14 Technical data). Ensure correct polarity ('+' and '-') when inserting batteries. Close the battery lid. Set date and time as described below.


Fig. 1

- If the Change battery symbol  is permanently illuminated, blood pressure can no longer be measured and you must replace all the batteries.

### Setting the date and time

Be sure to set the date and time correctly. This is the only way to save your measured values correctly with date and time for subsequent retrieval.

To switch to setting mode, reinsert the batteries or press and hold the SETTINGS button  for 5 seconds. Proceed as follows:

#### Measurement memory

The alarm flashes on the display.

- Set the alarm to off ('OFF') or on ('ON ) using the SETTINGS button . Confirm by pressing the SAVE button .

When the alarm is switched on, the alarm hour (a) and the alarm minutes (b) flash consecutively on the display.

- Depending on the display, use the SETTINGS button  to select the current alarm hour or alarm minutes and confirm using the SAVE button .
- When the alarm sounds, press any button to switch it off.

## Date

The display shows the year (c), the month (d) and the day (e).

- Depending on the display, use the SETTINGS button  to select the year, month or day and confirm using the SAVE button .

(c)

(d)

(e)

## Time

The hour (f) and the minutes (g) flash consecutively on the display.

- Depending on the display, use the SETTINGS button  to select the current current hour or minutes and confirm using the SAVE button .

(f)


(g)

To switch off the monitor, press the START/STOP button , otherwise the monitor will automatically switch off after 30 seconds.

## 6. Measuring blood pressure

### The 10 golden rules of blood pressure measurement

Many factors play a role in blood pressure measurement. These ten general rules help you to correctly carry out the measurement.


**1.** Rest for about 5 minutes before measuring your blood pressure. Even working at your desk increases blood pressure by an average of about 6 mmHg systolic and 5 mmHg diastolic.

**6.** Please do not talk or move during the measuring procedure. Talking increases values by about 6–7 mmHg.

**2.** Do not consume any nicotine or coffee for up to one hour before measurement.

**7.** Wait at least one minute between two measurements, so that the vessels are relieved from pressure in preparation for a new measurement.

**3.** Do not measure when you have a strong urge to urinate. A full bladder can increase blood pressure by about 10 mmHg.


**8.** The date and time are conveniently and easily recorded with the measured values as well as any medications you have taken with Veroyal® medi.connect.

**4.** Take measurements on the bare upper arm and while sitting upright.

**9.** Take measurements regularly. Even if your values have improved, you should continue to check them for monitoring purposes.

**5.** In the case of using a wrist monitor, hold the cuff at heart level during the measuring procedure. The cuff of an upper arm monitor is naturally positioned at the correct height.

**10.** Always take measurements at the same time of day. Because a person has about 100,000 different blood pressure values every day, individual measurements are not meaningful. Only regular measurements taken at the same time each day over a longer period of time give a meaningful assessment of blood pressure values.

## Applying the cuff

- Before applying the cuff, insert the cuff socket into the connector socket on the left side of the monitor.
- Do not mechanically constrict, compress or bend the sleeve tube.
- The blood pressure must be measured on a bare upper arm. If the cuff is completely open, thread the end of the cuff through the metal ring to form a loop. The hook-and-loop fastener must face outwards. Hold the cuff by the grip tab A (see Fig. 1) and pull it over the upper arm.


Fig. 1


Fig. 2

- The opening B (see Fig. 2) on the cuff, opposite the grip tab, should lie in the inner elbow. The tube should lie in the middle of the inner elbow and point towards the hand.
- Now bend your arm slightly, grasp the free end of the cuff, wrap it firmly round your arm and close the hook-and-loop fastener.

 The cuff should fit securely but not tightly. You should be able to insert two fingers between the arm and the cuff. Make sure the tube is not bent or damaged.


 Important: It is essential that the cuff is correctly applied to obtain correct measurements. The cuff included with the monitor is a universal cuff for upper arm circumferences of 22 cm to 42 cm. The white arrow must point to an area within the sizing scale. If it is outside the range, it is not possible to guarantee a correct measurement result.

 The innovative Veroval® monitor with Comfort Air technology ensures a comfortable measurement. Inflation to 190 mmHg is required for the first measurement. For subsequent measurements, inflation pressure is adjusted individually based on previously measured blood pressure values. This allows for comfortable upper arm measurements.

## Measuring the blood pressure

- You should take your blood pressure in a quiet place, in a relaxed and comfortable seated position.
- Measurement can be taken on the right or left arm. Over the long term, the arm giving higher readings should be used for blood pressure monitoring. However, if there is a very clear difference between readings on either arm, you need to check with your doctor which arm you should use for the measurement.
- Always measure on the same arm and place your forearm in a relaxed position on a support.
- We recommend that you measure your blood pressure while sitting with your back supported by the back of the chair. Place both feet flat on the floor next to each other. Do not cross your legs. Relax your forearm and hand with your palm facing upwards on a support and make sure the cuff is at the height of your heart.
- Do not take your blood pressure after taking a bath or playing sport.
- Do not eat, drink or exercise for at least 30 minutes before the measurement.
- Please wait for at least one minute between two measurements.
- The measured value is automatically stored in the displayed user memory. To change the user memory, press the SAVE button  when the device is switched off. To switch between user memories  and , hold the SAVE button  down for 3 seconds. Allocation to a different user memory is not possible after the measurement.
- Start a measurement only after applying the cuff. Press the START/STOP button  for about 2 seconds. The appearance of all display segments followed by the time indicates that the monitor is carrying out its self-test and is ready for use.

- Check the display segments for completion.


- After about 3 seconds the cuff is automatically inflated. If this inflation pressure is insufficient or if the measurement is interrupted, the device continues to pump in increments of 40 mmHg until a high enough pressure is reached.
- The increasing cuff pressure is displayed during inflation.

Important: You must not speak or move during the measurement.

- As the air pressure in the cuff decreases, the dropping pressure in the cuff is shown on the display.
- Once the measurement is complete, the systolic and diastolic blood pressure, pulse rate, memory position and results indicator are shown at the same time. After a few seconds, the date and time start to alternate with the memory position number. You can use the results indicator on the right side of the display to classify your measurement result (see Chapter 3 Blood pressure information). The measured value is automatically stored in the preset user memory.

- In addition to the measured values, the associated user memory  or  appear on the bottom right.
- To switch off the monitor, press the START/STOP button , otherwise the monitor will automatically switch off after 30 seconds.

If you wish to stop the measurement for any reason, simply press the START/STOP button.  The cuff inflation or measurement is interrupted and the pressure is automatically released.

### **How reliable is the blood pressure measurement with your device in the case of cardiac arrhythmia?**

Your blood pressure meter is generally suitable for measuring blood pressure even in the presence of cardiac arrhythmias, because the measurement algorithm filters these out.

In the case of severe cardiac arrhythmias, however, the blood pressure measurement may be affected because in the worst case there is permanent interference with the blood pressure signal.

If you occasionally suffer from cardiac arrhythmias, we recommend that you carry out an ECG diagnostic after each blood pressure measurement to completely rule out any measuring error. If the ECG confirms acute cardiac arrhythmia, repeat the blood pressure measurement after about 5–10 minutes.

## 7. Recording an ECG

There are two different methods used to record the ECG. Start with method 1 'Right index finger – chest'. If you cannot record an ECG, or the ECG is unstable (error message 'EE' frequently shown), change to method 2 'Left hand – right hand'.

Depending on the user's specific heart configuration (shape of the heart), the type or approach described in one of the methods indicated will be suitable. If you cannot record any stable measurements using a particular method, it may simply be due to harmless causes such as the heart shape but it may also be due to pathological reasons.

 Method 2 is very comfortable but the measurement stability is higher for method 1.

Ensure that there is no skin contact between your right and left hands (method 2) or your hand and chest (method 1). The recording will not be carried out correctly if this is the case. Do not move during the recording, do not speak and hold the device still. Movements of any type will distort the recordings.

Keep the pressure during the recording steady. Do not press the electrodes too firmly on the skin because the muscle tension that results may lead to imprecise measured values.

### **Method 1: Right index finger – chest**

- Hold the device in your right hand. Place your right index finger as illustrated on the upper electrodes (Fig. 1).

- Place the lower electrode of the device on bare skin about 5 cm below your left nipple. Hold the device pressed carefully against your skin until you hear a clicking sound (Fig. 2).


Fig. 1


Fig. 2

**Method 2: Left hand – right hand**

- Place your right index finger on the upper electrodes of the device. Place one finger from your left hand on the lower electrode (Fig. 3).
- Hold the device pressed carefully against your skin until you hear a clicking sound.



Fig. 3

**Carry out the recording**

- The device automatically switches on and **ECG** flashes briefly on the display.
- The device now starts to record and a 30-second countdown appears on the display. Continue to press the device against your skin until the 30 seconds has passed. During the recording your current heart rate will be shown on the display. A ♥ symbol also flashes in the same rhythm as your heartbeat.
- When the 30-second countdown has finished, the recording is complete and your final result appears on the display.

- To switch the device off after recording, press the START/STOP button **⏻** or wait 30 seconds until the device automatically switches off.


### ECG results display

The following results may be shown on the display after the recording.

	Results for a normal ECG recording.
	Indications of a slow heart rate (bradycardia) that is lower than 55 beats per minute.
	Indications of a fast heart rate (tachycardia) that is higher than 100 beats per minute.
	Indications of one or more pauses in the heart cycles that are longer than 2 seconds in each case.
	Indications of an arrhythmia during the ECG recording. Generally, the following arrhythmias are detected: Supraventricular arrhythmias (including atrial fibrillations / atrial flutters / pathological sinus arrhythmias / paroxysmal atrial tachycardia [if this persists, the message FAST appears] / supraventricular extrasystoles).
 	Indications of an arrhythmia during the ECG recording. Changed wave shape. Ventricular arrhythmias (including singular ventricular extrasystoles / bigeminy / trigeminy / series of ventricular extrasystoles / multifocal ventricular extrasystoles / ventricular tachycardia).
	Changed wave shape Ventricular arrhythmias

 It is recommended, particularly if status messages that deviate from  are repeatedly issued, to undergo a medical examination.

**i** If the displayed heart rate flashes, this means that the ECG signals are unstable or weak. In this case, repeat the recording.

**i** The contraction of the heart muscle is stimulated by electrical signals. A disorder of these electrical signals is called an arrhythmia. This may be due to a genetic predisposition, stress, age, lack of sleep, exhaustion or similar factors. A doctor can establish whether an irregular heartbeat is the result of an arrhythmia.

**Possible reasons for an invalid ECG recording**

The ECG recording can be invalid for the following reasons:

<p>The contact between the electrodes and the finger is not sufficiently firm.</p>	<p>Recordings that are made over or through clothing.</p>
<p>The device is held in the wrong direction and the wrong side is pressed against the chest.</p>	<p>The device is held in the left hand for method 1.</p>

## 8. Memory function

### User memory

- The Veroval® ECG and blood pressure monitor saves a total of up to 64 measurements in each user memory for ECG (Fig. 1) and/or blood pressure (Fig. 2), including the date and time of the measurement.
- Memory recall is activated by pressing the SAVE button  when the monitor is switched off. By repeatedly pressing the SAVE button , all stored values can be recalled one by one. To switch between user memories  and , hold the SAVE button  down for 3 seconds.


Fig. 1


Fig. 2

- You can cancel viewing the saved measured values at any time by pressing the START/STOP button . Otherwise, the device will switch off automatically after 30 seconds.
- Even if the power supply is interrupted, e.g. when changing the batteries, the stored values will still be available once full batteries are reinserted.

### Deleting saved readings

All data saved for each user can be separately deleted for user memory  and user memory . If you would like to delete the entire memory for a particular user, press the SAVE button  again and hold it down together with the START/STOP button  for 5 seconds. All data in the selected user memory are now deleted (Fig. 1). If you release the button ahead of time, no data will be deleted.


Fig. 1

## 9. Transferring measured values to Veroval® medi.connect

- Download the Veroval® medi.connect software from the website [www.veroval.info](http://www.veroval.info). Compatible with any computer with Windows 7, 8 or 10 (for the duration of official Microsoft support).

- Start the program and connect the Veroval® ECG and blood pressure monitor to your PC using the enclosed USB *Bluetooth*® adapter. To do this, insert the USB *Bluetooth*® adapter into the USB port of your computer. The Veroval® ECG and blood pressure monitor transfers the saved measured values via *Bluetooth*® directly to the Veroval® medi.connect software. Then follow the instructions in the Veroval® medi.connect software.

No data should be transferred during a measurement.

-  **Bluetooth**® is displayed on the ECG and blood pressure monitor.
- Start the data transfer using the 'medi.connect' software.
- During the data transfer, the data transfer symbol  flashes.

The ECG and blood pressure monitor switches off automatically after 30 seconds of inactivity and when communication with the PC is interrupted.

## 10. Error messages

Error that has occurred	Possible causes	Remedy
Monitor cannot be switched on	Batteries are missing, inserted incorrectly or flat.	Check batteries and insert four identical new batteries if necessary.
Cuff will not inflate	Cuff connector plug is incorrectly positioned in the connector socket of the monitor.	Check the connection between the cuff connector plug and socket.
	Wrong cuff type connected.	Verify that only an approved Veroval® cuff and connector are being used.
The pulse value flashes when the ECG measurement is displayed	Pulse could not be detected correctly.	Repeat the ECG measurement.
	Invalid ECG measurement	Check whether there is sufficient skin contact to the upper and lower electrodes when recording the ECG.

Error that has occurred	Possible causes	Remedy
E0	No pulse could be detected during the blood pressure measurement because the cuff was not correctly applied.	Check that the cuff has been positioned correctly. Apply the cuff so that two fingers fit between the cuff and the upper arm.
		Air tube is not correctly connected to the monitor. Check that the connector plug is correctly positioned. If this error occurs frequently, you should use a new cuff.
E1	The cuff pressure is greater than 300 mmHg.	Check that the cuff has been positioned correctly. Apply the cuff so that two fingers fit between the cuff and the upper arm.
	No pressure could be built up in the cuff.	Air tube is not correctly connected to the monitor. Check that the connector plug is correctly positioned. If this error occurs frequently, you should use a new cuff.
E2	Unusually high or low blood pressure measured values were detected.	Check that the cuff has been positioned correctly. Apply the cuff so that two fingers fit between the cuff and the upper arm.
		Air tube is not correctly connected to the monitor. Check that the connector plug is correctly positioned. If this error occurs frequently, you should use a new cuff.

Error that has occurred	Possible causes	Remedy
E3	The blood pressure measurement could not be carried out.	<p>Check that the cuff has been positioned correctly. Apply the cuff so that two fingers fit between the cuff and the upper arm.</p> <p>Air tube is not correctly connected to the monitor. Check that the connector plug is correctly positioned. If this error occurs frequently, you should use a new cuff.</p>
E4	The blood pressure measurement could not be carried out because the batteries are almost exhausted.	Replace batteries.
	The batteries are almost exhausted.	Replace batteries.
Implausible measured values during the blood pressure measurement	Implausible blood pressure values are often due to inappropriate handling of the monitor or mistakes being made when measuring the blood pressure.	Please observe the 10 golden rules of blood pressure measurement (see 53) and the safety information. Then repeat the measurement.

- Switch the monitor off if an error symbol appears. Check for all possible causes and note the 10 golden rules as well as the instructions for self-measurement in Chapter 2 Important information. Relax for at least 1 minute and then take the measurement again.
- The wordmark *Bluetooth*<sup>®</sup> and the accompanying logo are registered trademarks of Bluetooth SIG, Inc. Any use of these marks by PAUL HARTMANN AG is done under licence. Other trademarks and brand names are the property of the particular owner.

## 11. Caring for the device

- Only ever use a soft, damp cloth to clean the monitor. Please do not use thinners, alcohol, detergents or solvents. If the electrode surfaces of the device are soiled, clean them with a cotton bud soaked in cleaning alcohol.
- The cuff can be cleaned carefully using a slightly damp cloth and mild soap solution. Do not completely immerse the cuff in water.
- It is recommended to clean and disinfect the cuff regularly or after each use, especially when used by several users, to prevent infection. The cuff should be disinfected, particularly on the inside, by wiping with a disinfectant. Use a disinfectant that is compatible with the cuff materials. To protect from external influences, keep the monitor, the cuff and these instructions in the storage bag.

## 12. Warranty conditions

- We offer a 3-year warranty for this premium ECG and blood pressure monitor from the day of purchase and in accordance with the conditions listed below.
- Claims must be made during the warranty period. The date of purchase must be documented by an appropriately completed and stamped warranty certificate or proof of purchase.
- Within the warranty period, HARTMANN will replace or repair free of charge any faulty device components caused by material faults or manufacturing errors. This does not extend the warranty period.
- The device is intended for the purpose described in this instruction manual only.
- The warranty does not cover damage resulting from improper use or unauthorised interference. Accessory parts that are subject to wear and tear (batteries, cuffs, etc.) are excluded from the warranty. Claims for compensation are limited to the value of the merchandise; compensation for consequential damage or injury is expressly excluded.
- In warranty cases, please send the monitor with cuff and, if applicable, the power supply unit, together with the fully completed and stamped warranty certificate or purchase receipt directly or via your dealer to the Customer Services department for your country.

## 13. Contact information for customer enquiries

ZA HARTMANN South Africa  
Epsom Avenue, Northriding,  
2169 Johannesburg  
www.veroval.info  
Tel. +27 860 4278 6266  
phzahelpdesk@hartmann.info

Date of revision of the text: 2017-03

## 14. Technical data

Model:	Veroval® ECG and blood pressure monitor
Type:	BP750X
Measuring method:	Oscillometric blood pressure measurement, single-channel ECG ECG signal referenced to mass (earth)
Display range:	0–300 mmHg
Measuring range:	Systolic (SYS): 60–280 mmHg Diastolic (DIA): 30–200 mmHg Pulse: 30–180 pulse beats per minute We cannot guarantee that values outside the measuring range will be displayed correctly.
ECG bandwidth/sampling rate:	0.05 to 40Hz/256Hz
Display unit:	1 mmHg
Technical measurement accuracy:	Cuff pressure: $\pm 3$ mmHg Pulse: $\pm 5\%$ of the displayed pulse rate
Clinical measuring accuracy:	Complies with the requirements of DIN EN 1060-4 and DIN EN ISO 81060-2; Korotkoff validation method: Phase I (SYS), Phase V (DIA)
Operating mode:	Continuous mode
Power supply:	4× 1.5V alkaline manganese (AAA/LR03) bat- teries

Protection against electric shock:	Internal power supply ME device Applied part: Type BF
Protection against harmful ingress of water or solid materials:	IP22 (protected against penetration by solid objects with a diameter $\geq 12.5$ mm. Protected against water drops falling at an angle of up to $15^\circ$ to the monitor.)
Inflation pressure:	About 190 mmHg for first measurement
Automatic switch-off function:	30 seconds
Cuff:	Veroval® cuff for ECG and blood pressure monitor, secure fit cuff for arm circumference of 22 to 42 cm
Memory capacity:	2 × 64 measurements
Operating conditions:	Ambient temperature: $+10^\circ\text{C}$ to $+40^\circ\text{C}$ Relative humidity: 15–85%, non-condensing Air pressure 800–1050 hPa
Storage / transport conditions:	Ambient temperature: $-20^\circ\text{C}$ to $+55^\circ\text{C}$ Relative humidity: 15–85%, non-condensing Air pressure 800–1050 hPa
Serial number:	Inside the battery compartment
Computer interface to PC:	The USB <i>Bluetooth</i> ® adapter and the Veroval® medi.connect software enables the measured value memory to be read and a graphical representation of measured values to be displayed on a personal computer (PC).
Reference to standards:	IEC 60601-1; IEC 60601-1-2

**Technical data USB *Bluetooth*<sup>®</sup> adapter**

Interface Type:	Full-Speed USB Interface / with USB Type A Connector
Frequency Band:	2.4GHz
Networking <i>Bluetooth</i> <sup>®</sup> :	Bluetooth 4.0 Single Mode BLE
Antenna:	Internal
Display:	LED indicator
Operation Temperature/ Humidity:	+10°C to +40°C; 15~85%RH; 700-1060hPa
Storage Temperature/Humidity	-22°C to +55°C; 10~85%RH; 700-1060hPa
Power Source:	Notebook / PC
Dimensions:	L40,0 x W16,0 x H4,5 mm
Weight:	Approximately 2.4 g

**Legal requirements and guidelines**

- The Veroval<sup>®</sup> ECG and blood pressure monitor conforms to European regulations, which are based on the European Medical Device Directive 93/42/EEC, and bears the CE mark.
- The device complies, for example, with the European Standard EN 1060: Non-invasive sphygmomanometers - Part 1: General requirements and Part 3: Additional requirements for electro-mechanical blood pressure measuring systems and the standard EN 80601-2-30.
- The clinical testing of measuring accuracy was performed according to the European Standard EN 1060-4 as well as to EN 81060-2.
- Beyond statutory requirements, the monitor has been clinically validated by the ESH-IP2 Protocol of the European Society of Hypertension (ESH).
- The ECG device complies with the Medical Device Directive and the standards IEC 60601-2-25 (Medical electrical equipment – Part 2-25: Particular requirements for the safety of electrocardiographs) and IEC 60601-2-47 (Medical devices – Part 2-47: Particular requirements for the basic safety and essential performance of ambulatory electrocardiographic systems).
- We hereby guarantee that this product complies with the European R&TTE Directive 1999/5/EC.

## Electromagnetic compatibility

Table 1

**For all MEDICAL-ELECTRICAL (ME) DEVICES and ME SYSTEMS**

### Guidance and Manufacturer's Declaration – Electromagnetic emissions

The Veroval® ECG and blood pressure monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Veroval® ECG and blood pressure monitor should assure that it is used in such an environment.

<b>Emissions measurements</b>	<b>Concordance</b>	<b>Electromagnetic Environment – Guidance</b>
RF emissions CISPR 11	Group 1	The Veroval® ECG and blood pressure monitor uses RF energy exclusively for its internal functioning. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Veroval® ECG and blood pressure monitor is suitable for use in all establishments. The exception is domestic establishments and those directly connected to the public low-voltage network that supplies buildings used for domestic purposes.
Harmonic current emissions IEC 61000-3-2	Class A	
Emissions of voltage fluctuations/flicker	Complies	

**Table 2**

**For all MEDICAL-ELECTRICAL (ME) DEVICES and ME SYSTEMS**

**Guidance and Manufacturer’s Declaration – Electromagnetic immunity**

The Veroval® ECG and blood pressure monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Veroval® ECG and blood pressure monitor should assure that it is used in such an environment.

<b>Electromagnetic interference immunity tests</b>	<b>IEC 60601 test level</b>	<b>Compliance Level</b>	<b>Electromagnetic Environment – Guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If the floor has a synthetic material covering, the relative humidity must be at least 30%.
Rapid transient electrical interference / bursts in accordance with IEC 61000-4-4	± 2 kV for power supply cables	± 2 kV for power supply cables	The quality of the supply voltage should be that of a typical business or hospital environment.
Surge voltages (surges) in accordance with IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth(s)	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth(s)	The quality of the supply voltage should be that of a typical business or hospital environment.
Voltage drops, brief interruptions and fluctuations in the supply voltage in accordance with IEC 61000-4-11	< 5% $U_T$ (> 95% drop in $U_T$ ) for 0.5 period	< 5% $U_T$ (> 95% drop in $U_T$ ) for 0.5 period	The quality of the supply voltage should be that of a typical business or hospital environment. If the user of the Veroval® ECG and blood pressure monitor requires continued operation during a power mains interruption, it is recommended that the Veroval® ECG and blood pressure monitor be powered from an uninterruptible power supply or battery.
	40% $U_T$ (60% drop in $U_T$ ) for 5 periods	40% $U_T$ (60% drop in $U_T$ ) for 5 periods	
	70% $U_T$ (30% drop in $U_T$ ) for 25 periods	70% $U_T$ (30% drop in $U_T$ ) for 25 periods	
	< 5% $U_T$ (> 95% drop in $U_T$ ) for 5 sec.	< 5% $U_T$ (> 95% drop in $U_T$ ) for 5 sec.	
Magnetic field at supply frequency (50/60 Hz)	3 A/m	Not applicable	Not applicable

**NOTE:**  $U_T$  is the mains alternating voltage before use of the test levels.

**Table 3**

**For all NON-LIFE-SUPPORT MEDICAL-ELECTRICAL (ME) DEVICES and ME SYSTEMS. Guidance and Manufacturer’s Declaration – Electromagnetic immunity**

The Veroval® ECG and blood pressure monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Veroval® ECG and blood pressure monitor should assure that it is used in such an environment.

Electromagnetic interference immunity tests	IEC 60601 test level	Compliance Level	Electromagnetic Environment – Guidance
<p>Conducted RF disturbances IEC 61000-4-6</p> <p>Radiated RF disturbances IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile radio devices should not be used closer to the Veroval® ECG and blood pressure monitor, including the cables, than the recommended safety distance, which is calculated according to the equation applicable to the transmission frequency.</p> <p>Recommended separation distance</p> $d = 1.2 \sqrt{P}$ <p><math>d = 1.2 \sqrt{P}</math> for 80 MHz to 800 MHz</p> <p><math>d = 2.3 \sqrt{P}</math> for 800 MHz to 2.5 GHz</p> <p>Where <math>P</math> is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).</p>
<p><b>NOTE 1:</b> At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p><b>NOTE 2:</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Veroval® ECG and blood pressure monitor is used exceeds the applicable RF compliance level above, the Veroval® ECG and blood pressure monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Veroval® ECG and blood pressure monitor.</p>			

**Table 4**

**For all NON-LIFE-SUPPORT MEDICAL-ELECTRICAL (ME) DEVICES and ME SYSTEMS. Recommended separation distances between portable and mobile RF communications equipment and the Veroval® ECG and blood pressure monitor.**

The Veroval® ECG and blood pressure monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the Veroval® ECG and blood pressure monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Veroval® ECG and blood pressure monitor, as recommended below, according to the maximum output power of the communications equipment.

Nominal output of the transmitter W	Protective distance dependent on the transmitter frequency m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Made under the control of Paul Hartmann AG by
MD Biomedical, Inc.
Medical Device Safety Service GmbH